

## Pulmonary CTA enhanced with blood flow maps. Where there is smoke there is fire.

*"<sup>SURE</sup>Subtraction Lung enhances routine pulmonary CTA studies with color iodine maps for the evaluation of pulmonary perfusion. Overall diagnostic confidence is greatly improved particularly to identify perfusion deficit in the case of small or occult pulmonary emboli."*

**Dr. Jean-François Paul**


Chief of Cardiac Imaging, Department of Radiology  
Institut Mutualiste Montsouris, Paris, France



### Patient history

A 63-year-old female presented for suspected Pulmonary Embolism. The scan was performed with 35 ml of contrast and iodine maps were created automatically.

### Results




The color iodine maps show a large wedge-shaped perfusion defect in the middle lobe of the right lung. Further inspection revealed an embolic occlusion of a distal branch of the right pulmonary artery.

## Technology

Utilizing a pre- and post-contrast scan, the <sup>SURE</sup>Subtraction Lung application isolates the iodine signal and displays the result as a color overlay. A dedicated scan mode synchronizes the pre- and post-contrast scans and automatically outputs iodine maps directly to the reading station with no need for additional user interaction.

The key to obtaining accurate results lies with an anatomically aware 3D deformable registration algorithm that compensates for patient motion that may occur between two scans.



## Conclusion

The addition of iodine maps to pulmonary CTA enables the assessment of the distribution of contrast media in the pulmonary parenchyma. Subtraction CT of the lungs can automatically provide iodine maps in every routine pulmonary CTA examination. Similar results were reported by Grob et al.\*

\* Grob, D., Oostveen, L.J., Prokop, M. et al. Imaging of pulmonary perfusion using subtraction CT angiography is feasible in clinical practice. Eur Radiol (2018).  
<https://doi.org/10.1007/s00330-018-5740-4>

## CANON MEDICAL SYSTEMS CORPORATION

<https://global.medical.canon>

©Canon Medical Systems Corporation 2019. All rights reserved.  
Design and specifications are subject to change without notice.  
Model number: TSX-305A MCACT0335EA 2019-03 CMSC/Produced in Japan

Canon Medical Systems Corporation meets internationally recognized standards for Quality Management System ISO 9001, ISO 13485. Canon Medical Systems Corporation meets the Environmental Management System standard ISO 14001.

Aquilion ONE, Aquilion ONE GENESIS, <sup>SURE</sup>Subtraction, <sup>SURE</sup>Exposure and Made for Life are trademarks of Canon Medical Systems Corporation.

Clinical results may vary due to clinical setting, patient presentation and other factors.

## Acquisition

|                 |  |
|-----------------|--|
| Scanner Model:  | Aquilion ONE / GENESIS Edition                 |
| Scan Mode:      | Ultra Helical                                  |
| Collimation:    | 0.5 mm × 80                                    |
| Exposure:       | 100 kV<br><sup>SURE</sup> Exposure             |
| Rotation Time:  | 0.275 second                                   |
| Dose Reduction: | AIDR*1 3D Enhanced                             |
| CTDI:           | 0.7 mGy pre-contrast/<br>3.3 mGy post-contrast |
| DLP:            | Total 155.1 mGy-cm                             |
| Effective Dose: | 2.17 mSv                                       |
| k-factor:       | 0.014*2  |
| Application:    | <sup>SURE</sup> Subtraction Lung*3             |

\*1 Adaptive Iterative Dose Reduction

\*2 American Association of Physicists in Medicine (AAPM) Report 96, 2008.

\*3 Option

*Made For life*