

Prospective Clinical Evaluation of Novel Deep Learning Denoising Reconstruction for Improved Image Quality in Musculoskeletal Magnetic Resonance Imaging

Hung P. Do, PhD

Manager Medical Affairs – Sr. Clinical Scientist
Canon Medical Systems USA, Inc.

Dawn Berkeley

Manager Medical Affairs – Sr. Clinical Development
Canon Medical Systems USA, Inc.

Introduction

Over the last few years, deep learning has been widely applied to magnetic resonance image reconstruction (1-10). It has been demonstrated that deep learning-based reconstruction (DLR) may provide performance gains when compared to conventional reconstruction methods (11). In MRI, a high resolution and signal-to-noise ratio (SNR) are desired clinically for better delineation of subtle anatomical and pathological features that are critical in making a clinical diagnosis. Higher resolution and SNR require longer scan times, which is often less desired in a clinical environment. Deep learning techniques may have great potential. Through the denoising process, image quality and SNR gain could assist clinicians in providing flexibility in scan parameters and be used to improve resolution and/or shorten scan time.

In this paper, Advanced Clear-IQ intelligent Engine (AiCE) DLR was prospectively evaluated, in which the results suggest that AiCE DLR gains in image quality, denoising, and SNR may help clinicians in providing flexibility in scan parameters and be used to improve resolution and/or shorten scan time in comparison to images routinely acquired at the Institution.

Evaluation in a Clinical Setting

This prospective clinical study was performed in a clinical site in the USA, in which images acquired from 137 sequences

in 40 patients (10 hips, 10 shoulders, 20 knees) were collected. The study was carried out on a 3T MRI scanner (Canon Vantage Galan) with approved IRB. Each sequence was imaged twice: first with a standard clinical protocol (reference sequence, or “REF”), then with a modified protocol that improved resolution by 48% and shortened scan time by 8%. The modified protocol was reconstructed with Advanced intelligent Clear-IQ Engine (AiCE) DLR (“DLR”) (10) and three other clinically used reconstruction filters (named NL2, GA43, and GA53) resulting in a total of five reconstructions for each sequence. As a blinded, randomized clinical image review, all reconstructions were de-identified and randomized before sharing with three experienced American Board-Certified MSK radiologists for review. Representative images for the reference images and images reconstructed with different reconstruction methods are shown in Figure 1, in which arrows indicate pathologies identified by one of the radiologists.

For image quality assessment, radiologists were asked to assign scores of image quality based on the entire image and on specific anatomical and pathological features they observe in the image (Table 1). The reviewers were also asked to rank the images in order of preference for different reconstruction methods, in which 5 means the most preferred and 1 means the least preferred.

The average reader score was highest for DLR across all eight scoring criteria and the forced ranking compared to the other reconstruction methods (Figure 2).




Figure 1: Representative images acquired with the standard clinical protocol (reference image, "REF") and those acquired using a modified protocol for higher resolution and/or shorter scan time and reconstructed with DLR and conventional reconstructions (NL2, GA43, and GA53). Arrows represent pathologies identified by one of the radiologists.

| Score | Overall Image Quality Assessment | | | | Assessment of Specific Anatomical & Pathological Features | | | |
|-------|---|---|---|------------------------------------|---|--|-----------------------------------|--|
| | Overall Image Noise | Overall Image Sharpness | Overall Image Degradation/ Artifacts | Overall Diagnostic Confidence | Features Sharpness | Features Contrast | Features Conspicuity | Features Diagnostic Confidence |
| 1 | All structures are too noisy | All structures are not sharp on most images | All structures are degraded | Diagnostic confidence is very poor | Sharpness of features is very poor | Contrast between features and surrounding is very poor | Features conspicuity is very poor | Diagnostic confidence for the reviewed features is very poor |
| 2 | Most structures are too noisy | Most structures are not sharp on some images | Most structures are degraded | Poor | Poor | Poor | Poor | Poor |
| 3 | A few structures are noisy on most images | Most structures are sharp on some images | A few structures are degraded on most images | Fair | Fair | Fair | Fair | Fair |
| 4 | A few structures are noisy on a few images | Most structures are sharp on most images | A few structures are degraded on a few images | Good | Good | Good | Good | Good |
| 5 | There is no appreciable noise on any of the relevant images | All structures are sharp on all relevant images | There is no appreciable image degradation on any of the relevant images | Excellent | Excellent | Excellent | Excellent | Excellent |

Table 1: Scoring criteria and scoring instructions




Figure 2: Average scores of three musculoskeletal radiologist readers for the five reconstruction methods. DLR was rated significantly higher across all eight scoring criteria and the forced ranking compared to the other reconstruction methods.

Summary

AiCE DLR removes noise, hence improved SNR.* The SNR improvement may help clinicians in providing flexibility in scan parameters and be used to help alleviate the fundamental tradeoffs between SNR, resolution, and scan time. This prospective clinical evaluation showed that AiCE DLR may enable clinicians to increase resolution and/or shorten scan time while allowing preferred image quality.

References

1. Yang Y, Sun J, Li H, Xu Z. Deep ADMM-Net for Compressive Sensing MRI. *Adv Neural Inf Process Syst*. Curran Associates, Inc.; 2016. p. 10–18.
2. Lee D, Yoo J, Ye JC. Deep residual learning for compressed sensing MRI. 2017 IEEE 14th Int Symp Biomed Imaging ISBI 2017. 2017. p. 15–18. doi: 10.1109/ISBI.2017.7950457.
3. Hammernik K, Klatzer T, Kobler E, et al. Learning a variational network for reconstruction of accelerated MRI data. *Magn Reson Med*. 2018;79(6):3055–3071. doi: 10.1002/mrm.26977.
4. Zhu B, Liu JZ, Cauley SF, Rosen BR, Rosen MS. Image reconstruction by domain-transform manifold learning. *Nature*. 2018;555(7697):487–492. doi: 10.1038/nature25988.
5. Akçakaya M, Moeller S, Weingärtner S, Uğurbil K. Scan-specific robust artificial-neural-networks for k-space interpolation (RAKI) reconstruction: Database-free deep learning for fast imaging. *Magn Reson Med*. 2019;81(1):439–453. doi: 10.1002/mrm.27420.
6. Chen F, Taviani V, Malkiel I, et al. Variable-Density Single-Shot Fast Spin-Echo MRI with Deep Learning Reconstruction by Using Variational Networks. *Radiology*. Radiological Society of North America; 2018;289(2):366–373. doi: 10.1148/radiol.2018180445.
7. Knoll F, Murrell T, Sriram A, et al. Advancing machine learning for MR image reconstruction with an open competition: Overview of the 2019 fastMRI challenge. *Magn Reson Med*. 2020;84(6):3054–3070. doi: 10.1002/mrm.28338.
8. Yaman B, Hosseini SAH, Moeller S, Ellermann J, Uğurbil K, Akçakaya M. Self-supervised learning of physics-guided reconstruction neural networks without fully sampled reference data. *Magn Reson Med*. 2020;84(6):3172–3191. doi: 10.1002/mrm.28378.
9. Akçakaya M, Yaman B, Chung H, Ye JC. Unsupervised Deep Learning Methods for Biological Image Reconstruction. *ArXiv210508040 Phys*. 2021; <http://arxiv.org/abs/2105.08040>.
10. Kidoh M, Shinoda K, Kitajima M, et al. Deep Learning Based Noise Reduction for Brain MR Imaging: Tests on Phantoms and Healthy Volunteers. *Magn Reson Med Sci*. 2020;19(3):195–206. doi: 10.2463/mrms.mp.2019-0018.
11. Knoll F, Murrell T, Sriram A, et al. Advancing machine learning for MR image reconstruction with an open competition: Overview of the 2019 fastMRI challenge. *Magnetic Resonance in Medicine*. 2020;84(6):3054–3070. doi:10.1002/mrm.28338.

The clinical results described in this paper are the experience of the authors. Results may vary due to clinical setting, patient presentation and other factors.

* AiCE provides higher SNR compared to typical low pass filters.

CANON MEDICAL SYSTEMS USA, INC.

<https://us.medical.canon>

2441 Michelle Drive, Tustin, CA 92780 | 800.421.1968

©Canon Medical Systems, USA 2023. All rights reserved.
Design and specifications are subject to change without notice.

Made for Life is a trademark of Canon Medical Systems Corporation.

MRWP14285US

Made For life