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2B702-067EN*T

CANON MEDICAL SYSTEMS CORPORATION

© CANON MEDICAL SYSTEMS CORPORATION 2015-2021. ALL RIGHTS RESERVED.

Guidelines for cleaning, disinfection, and sterilization of transducers

Issued: 2021-04

This manual describes the cleaning, disinfection, and sterilization procedures for the ultrasound transducer.
For the operating precautions and procedures for the transducers, refer to the operation manual for each transducer.

Trademarks

This manual may include trademarks or registered trademarks of companies other than Canon Medical Systems.

1. Safety Precautions

1.1 Meaning of Signal Words

In this manual, the signal words **DANGER**, **WARNING**, **CAUTION**, and *NOTICE* are used regarding safety and other important instructions. The signal words and their meanings are defined as follows. Please understand their meanings clearly before reading this manual.

Signal word	Meaning
⚠DANGER	Indicates an imminently hazardous situation which, if not avoided, will result in death or serious injury.
⚠WARNING	Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.
⚠CAUTION	Indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury.
<i>NOTICE</i>	Indicates a potentially hazardous situation which, if not avoided, may result in property damage.

1.2 Safety Precautions

Observe the following precautions to ensure the safety of patients as well as operators when performing cleaning, disinfection, or sterilization of the transducer.

⚠WARNING Never immerse the transducer connector or any other non-waterproof sections into liquids such as water or cleaning solution. Immersion may cause electric shock. Refer to the tables and figures in section 3 for the immersible range of each transducer model.

⚠CAUTION

1. Precautions concerning cleaning, disinfection, and sterilization.
 - 1) Observe the following precautions to prevent infection.
 - Wear protective gloves when performing cleaning.*
 - Wear sterile protective gloves when performing disinfection or sterilization.*
 - Wear new protective gloves each time cleaning, disinfection, and sterilization are performed.
 - Clean the transducer before and after examination. Disinfect or sterilize the transducer as required.
 - Note that local regulations may require that the transducer be disinfected before sterilization.
 - Sterilize the transducer and biopsy adaptor before and after an ultrasound-guided biopsy procedure is performed. Failure to do so may result in the transducer and biopsy adaptor becoming sources of infection.

* Refer to the FDA's March 29, 1991 Medical Alert on Latex Products.

⚠ CAUTION

- 2) After cleaning, rinse the transducer thoroughly with purified water to remove all chemical residues. After disinfection, rinse the transducer thoroughly with sterile or deionized water to remove all chemical residues.
Chemical residues on the transducer may be harmful to the human body.
- 3) After chemical cleaning or chemical disinfection, thoroughly dry the transducer surface.
- 4) After sterilization, degas the transducer. Gas residues on the transducer may be harmful to the human body.
- 5) The efficacy of the cleaning solutions, disinfectants, and sterilizing gases is not guaranteed by Canon Medical Systems. Contact the manufacturers for information on the activity of the products.
- 6) To ensure the prevention of infection, confirm the effectiveness of each chemical for cleaning, disinfection, or sterilization based on the criteria (such as effective period, number of times of use, discoloration, and results of using the effectiveness test kit) described in the documentation provided by the relevant manufacturer.
- 7) Do not clean, disinfect, or sterilize the transducer using chemicals or methods other than those specified in this guide. If chemicals or methods other than those specified in this guide are used, the transducer may not be properly cleaned, disinfected, or sterilized or may be damaged.
- 8) Observe the following precautions to prevent transducer malfunction.
 - Only soft materials such as soft cloth or soft gauze should be used when removing ultrasound gel from the transducer or wiping water or disinfectant from the transducer surface after cleaning or disinfection. Use of hard or abrasive cloth or gauze may damage the transducer.
 - The transducer must not be immersed in a chemical solution for more than three hours.
 - Do not permit the transducer to become overheated (more than 60°C (140°F)) during cleaning, disinfection, and sterilization.
 - The cleaning, disinfection, or sterilization conditions, such as the temperature and pressure, differ depending on the product. In addition, some products cannot be subjected to disinfection or sterilization procedures. Confirm the detailed conditions by referring to the tables and figures in section 3.

NOTE The mouthpiece supplied with the TEE transducer can be high-pressure steam sterilized (temperature: 134°C (273.2°F), holding time: 18 min). Confirm that "AUTOCLAVABLE" is indicated on the mouthpiece.

1.3 Chemical Hazard (for USA only)

Observe the following instruction in order to protect patients and operators from inflammation or poisoning by chemical substances.

⚠️ WARNING: This product can expose you to chemicals including lead and lead compounds, which are known to the State of California to cause cancer, and phthalates, which are known to the State of California to cause birth defects or other reproductive harm.
For more information go to www.p65warnings.ca.gov.

1.4 Precautions Concerning Restriction of Sale or Use (for USA only)

United States law restricts this device to sale or use by, or on the order of a physician.

2. Cleaning, Disinfection, Sterilization

This section describes the methods for cleaning, disinfection, and sterilization.

Cleaning and disinfection between patients must be performed to prevent transmission of disease.

All transducers must be thoroughly cleaned prior to disinfection.

The level of disinfection required is based on patient contact.

2.1 Cleaning

<<Items to be used: Protective gloves, cleaning solution or cleaning wipes, purified water, clean soft cloth or gauze, single-use sponge*>>

* The single-use sponge must not include any abrasive parts or contain any abrasive cleanser.

- (1) Wear protective gloves to prevent infection.
Wear new protective gloves each time cleaning is performed.
- (2) If an accessory that can be disassembled (e.g., biopsy adaptor) is provided, disassemble it. For details concerning such parts, refer to the operation manual for the transducer.
- (3) Wash off all organic materials (such as blood or other bodily fluids) from the transducer under purified water. A single-use sponge can be used for washing. Do not use a brush, because it may damage the transducer.
- (4) In accordance with the tables and figures in section 3, immerse the transducer in a cleaning solution or wipe the transducer using wipes to dissolve or remove all remaining organic materials. Use a single-use sponge if necessary. If dried organic materials are present on the transducer, immerse it in the cleaning solution for a prolonged period.
- (5) Remove all residual organic materials and cleaning solution from the transducer by rinsing it under purified water. Repeat steps (3) and (4) to visually confirm that all residual organic materials and cleaning solution have been completely removed from the transducer. Do not reuse the purified water.
- (6) Dry the surface of the transducer using clean soft cloth or gauze. Do not use heat to dry the transducer.
Handle the cleaning solution or wipes as described in the documentation provided by the relevant manufacturer. To maintain the effectiveness of the cleaning solution or wipes, ensure that the concentration, temperature, and other conditions specified in the documentation provided by the manufacturer are met. To confirm the effectiveness of the cleaning solution or wipes, use the criteria (such as effective period, number of times of use, discoloration, and results of using the effectiveness test kit) described in the documentation provided by the manufacturer.
- (7) Confirm that the transducer shows no signs of damage, deformation, or peeling.

2.2 Disinfection

Before disinfection, the transducer must be cleaned.

Note that local regulations may require that the transducer be disinfected before sterilization.

<<Items to be used: Sterile protective gloves, disinfectant, sterile water or deionized water, sterile soft cloth or gauze>>

- (1) Wear sterile protective gloves to prevent infection.
Wear new sterile protective gloves each time disinfection is performed.
- (2) Disinfect the transducer using the chemicals listed in section 3 "List of Chemicals". Ensure that the entire immersible range of the transducer is completely in contact with the disinfectant.
- (3) Rinse the transducer thoroughly with sterile or deionized water, as described in the documentation provided by the manufacturer, in order to remove all residual disinfectant. (Rinsing is not necessary when gas disinfection is performed using Trophon EPR or Trophon 2.) Do not reuse the sterile or deionized water.
- (4) Dry the surface of the transducer using sterile soft cloth or gauze.
Do not use heat to dry the transducer.
Handle the disinfectant as described in the documentation provided by the relevant manufacturer.
To maintain the effectiveness of the disinfectant, ensure that the concentration, temperature, and other conditions specified in the documentation provided by the manufacturer are met. To confirm the effectiveness of the disinfectant, use the criteria (such as effective period, number of times of use, discoloration, and results of using the effectiveness test kit) described in the documentation provided by the manufacturer.
- (5) Confirm that the transducer shows no signs of damage, deformation, or peeling.

2.3 Sterilization

Before sterilization, the transducer must be cleaned.

Note that local regulations may require that the transducer be disinfected before sterilization.

<<Items to be used: Sterile protective gloves, sterilant>>

* Some types of transducers cannot be sterilized or the sterilization conditions may differ.

- (1) Wear sterile protective gloves to prevent infection.
Wear new sterile protective gloves each time sterilization is performed.
- (2) Sterilize the transducer using the chemicals listed in section 3 "List of Chemicals".
* Place the transducer in a sterilization packing case and then place it in the sterilizer.
- (3) After gas sterilization, perform aeration to remove all gas residues on the transducer surface.
- (4) Confirm that the transducer shows no signs of damage, deformation, or peeling.


3. List of Chemicals

This section lists the chemicals used when washing, disinfecting, or sterilizing the transducer. Refer to the lists for the available chemicals.


In addition, refer to the manuals for the chemicals for the handling details.

- NOTE
1. When multiple model names included in the List of Chemicals are referred to, note that they are abbreviated. For example, "PLT-604AT and PLT-704AT" is abbreviated as "PLT-604AT/704AT".
 2. For the transducers supported by each diagnostic ultrasound system, refer to the operation manual supplied with the diagnostic ultrasound system.


- Type A (PC-20M etc.)


- Type B (PET-508MA etc.)

- Type C (PET-511BTM etc.)

- Type D (PLT-704AT etc.)

Never immerse the non-waterproof sections of the transducer into liquids such as water or cleaning solution. Immersion may cause electric shock.

Cleaning (1/2)

Country of origin :
 Any : Any country
 FRA : France
 DEU : Germany
 GBR : United Kingdom
 USA : United States
 JPN : Japan
 AUS : Australia
 CAN : Canada

Chemical name/type	Enzyme				
	Trade name	CIDEZYME®	3M™ Rapid Multi-Enzyme cleaner	ANIOSYME DD1 *2	neodisher® mediclean *2
Manufacturer	J&J	3M Company	Laboratories ANIOS	DR. WEIGERT GmbH & Co. KG	
Country of origin	Any	Any	FRA	DEU	
Concentration (Dilution ratio)	0.8% (125)	1% (100)	0.5% (200)	1% (100)	
Time	1 min.	10 min.	5 min.	10 min.	
Temperature	Room temperature		Room temperature		
Humidity	Normal humidity		Normal humidity		
Pressure	Normal pressure		Normal pressure		
Model name	Immersible range				
PC-20M	Type A	OK	OK	OK	OK
PC-50M	Type A	OK	OK	OK	OK
PET-508MA	Type B	OK	OK	-	-
PET-512MA/512MC	Type B	OK	OK	-	-
PET-511BTM*1	Type C	OK	OK	OK	-
PET-512MD*1, PEI-512VX*1	Type C	OK	OK*2	OK	-
PET-609MA*1	Type C	OK	OK*2	OK	-
PET-805LA	Type C	OK	-	-	-
PLT-604AT/704AT/805AT/1202S/1204BT PVT-375BT/375SC/382BT/661VT/674BT	Type D	OK	OK	OK	OK
PST-25BT/65AT	Type D	OK	OK	OK	OK
PVT-745BTF/745BTH	Type D	OK	OK	OK	OK
PLT-308P/308BTP, PVT-770RT	Type D	OK	OK	OK	OK
PLT-705BTF/705BTH	Type D	OK	OK	OK	OK
PVT-745BTV	Type D	OK	OK	OK	OK
PVT-350BTP	Type D	OK	OK	OK	OK
PLT-1204BX	Type D	OK	OK	OK	OK
PLT-704SBT/705BT/1005BT PVT-781VT/781VTE	Type D	OK	OK	OK	OK
PLT-704ST	Type D	OK	OK	OK	OK
PLT-1204ST, PVT-375ST	Type D	OK	OK	OK	OK

OK : Use of the chemical is permitted. : Use of the chemical is not permitted.

*1: If a waterproof cover is placed over the transducer connector, the entire transducer can be immersed in water or cleaning solution.

*2: Not for use in the USA.

Enzyme					Isopropyl alcohol/ Ethylene glycol monobutyl ether	Alkyldimethylbenzyl- ammonium chloride/2-ethanol/ Alkyl polysaccharide/ Ethylenediaminetetraacetic acid
Tristel Pre-Clean Wipes *2	MetriZyme™	Klenzyme®	Instru-Zyme® Gel	Endozime® Xtreme Power *2	CaviWipes™	Intercept® Wipes
Tristel Solutions Limited	Metrex Research	STERIS Corporation	Summit Medical	RUHOF Corporation	Metrex Research	MEDIVATORS Inc.
GBR	USA	USA	USA	USA	USA	USA
Working solution	0.77% (130)	0.8% (125)	Working solution	0.2% (500)	Working solution	Working solution
Wiping only	5 min.	5 min.	Wiping only	2 min.	Wiping only	
Room temperature		Room temperature			Room temperature	
Normal humidity		Normal humidity			Normal humidity	
Normal pressure		Normal pressure			Normal pressure	
OK	OK	OK	OK	OK	OK	OK
OK	OK	OK	OK	OK	-	-
OK	OK	-	-	OK	-	-
OK	OK	-	-	OK	-	-
OK	OK	OK	-	OK	-	-
OK	OK*2	-	-	OK	-	-
OK	OK*2	-	-	OK	-	-
-	-	-	-	OK	-	-
OK	OK	OK	OK	OK	OK	OK
OK	OK	OK	OK	OK	OK	OK
OK	OK	OK	OK	OK	OK	OK
OK	OK	OK	OK	OK	OK	OK
OK	OK	OK	OK	OK	OK	OK
OK	OK	OK	OK	OK	OK	OK
OK	OK	OK	OK	OK	OK	OK
-	OK	-	OK	OK	-	-
OK	OK	OK	OK	OK	OK	OK
OK	OK	OK	OK	OK	OK	OK
OK	OK	OK	OK	OK	OK	OK
OK	OK	OK	OK	OK	OK	OK
OK	OK	OK	OK	OK	OK	OK

Cleaning (2/2)

Country of origin :
 Any : Any country
 FRA : France
 DEU : Germany
 GBR : United Kingdom
 USA : United States
 JPN : Japan
 AUS : Australia
 CAN : Canada

Chemical name/type	Enzyme				
	Trade name	CIDEZYME®	3M™ Rapid Multi-Enzyme cleaner	ANIOSYME DD1 *2	neodisher® mediclean *2
Manufacturer	J&J	3M Company	Laboratories ANIOS	DR. WEIGERT GmbH & Co. KG	
Country of origin	Any	Any	FRA	DEU	
Concentration (Dilution ratio)	0.8% (125)	1% (100)	0.5% (200)	1% (100)	
Time	1 min.	10 min.	5 min.	10 min.	
Temperature	Room temperature		Room temperature		
Humidity	Normal humidity		Normal humidity		
Pressure	Normal pressure		Normal pressure		
Model name	Immersible range				
PST-25SX	Type D	OK	OK	–	OK
PST-25ST/50BT, PVT-712BT	Type D	OK	OK	OK	OK
PST-30BT	Type D	OK	OK	OK	OK
PLT-1204MV PVT-375MV/382MV/675MV/675MVL	Type D	OK	OK	OK	OK
PVT-681MV	Type D	OK	–	–	–
PVT-681MVL	Type D	OK	OK	OK	OK
PVL-715RST	Type D	OK	OK	OK	OK
PVI-475BX/475BT/574BX/482BX/450BXP PVT-475BT/482BT/574BT PLI-1205BX/2004BX/705BX/605BX PSI-30BX, PST-28BT/65BT	Type D	OK	OK	OK	OK
PSI-70BT	Type D	OK	OK	OK	OK
PSI-30VX/40VX	Type D	OK	OK	OK	OK
PSI-50VX	Type D	OK	OK	OK	OK
PLI-2002BT	Type D	OK	OK	OK	OK
PLI-3003BX	Type D	OK	OK	OK	OK
PLT-1202BT	Type D	OK	OK	OK	OK
PVT-675MVS	Type D	OK	OK	OK	OK

OK : Use of the chemical is permitted. : Use of the chemical is not permitted.

*1: If a waterproof cover is placed over the transducer connector, the entire transducer can be immersed in water or cleaning solution.

*2: Not for use in the USA.

Enzyme					Isopropyl alcohol/ Ethylene glycol monobutyl ether	Alkyldimethylbenzyl- ammonium chloride/2-ethanol/ Alkyl polysaccharide/ Ethylenediaminetetraacetic acid
Tristel Pre-Clean Wipes *2	MetriZyme™	Klenzyme®	Instru-Zyme® Gel	Endozime® Xtreme Power *2	CaviWipes™	Intercept® Wipes
Tristel Solutions Limited	Metrex Research	STERIS Corporation	Summit Medical	RUHOF Corporation	Metrex Research	MEDIVATORS Inc.
GBR	USA	USA	USA	USA	USA	USA
Working solution	0.77% (130)	0.8% (125)	Working solution	0.2% (500)	Working solution	Working solution
Wiping only	5 min.	5 min.	Wiping only	2 min.	Wiping only	
Room temperature		Room temperature			Room temperature	
Normal humidity		Normal humidity			Normal humidity	
Normal pressure		Normal pressure			Normal pressure	
-	OK	-	OK	OK	-	-
OK	OK	OK	OK	OK	OK	OK
OK	OK	OK	OK	OK	OK	OK
-	OK	-	OK	OK	OK	OK
-	-	-	-	OK	OK	OK
OK	OK	-	OK	OK	OK	OK
OK	OK	OK	OK	OK	OK	OK
OK	OK	OK	-	OK	-	-
OK	OK	OK	OK	OK	OK	OK
OK	OK	OK	-	OK	OK	-
OK	OK	OK	-	OK	OK	-
OK	OK	OK	-	OK	-	-
OK	OK	OK	-	OK	OK	-
OK	OK	OK	-	OK	OK	-

Low/Middle-Level Disinfection -1 (1/2) (Except for USA)

Note that local regulations may require that the transducer be disinfected before sterilization.

Country of origin :
 Any : Any country
 FRA : France
 DEU : Germany
 GBR : United Kingdom
 USA : United States
 JPN : Japan
 AUS : Australia
 CAN : Canada

Chemical name/type	Didecyldimethyl-ammonium chlorides/Quaternary ammonium compounds	Didecyldimethyl-ammonium chlorides	Glyoxal/Glutaral/Didecyldimethyl-ammonium chlorides
Trade name	CLEANISEPT® WIPES	WIP'ANIOS Excel	Lysoformin® 3000
Manufacturer	Dr. Schumacher GmbH	Laboratories ANIOS	Lysoform Disinfection Ltd.
Country of origin	DEU	FRA	DEU
Concentration (Dilution ratio)	Working solution	Working solution	2% (50)
Time	Wiping only	Wiping only	15 min.
Temperature	Room temperature	Room temperature	Room temperature
Humidity	Normal humidity	Normal humidity	Normal humidity
Pressure	Normal pressure	Normal pressure	Normal pressure
Model name	Pressure	Normal pressure	Normal pressure
PC-20M	Immersible range	OK	OK
	Type A	OK	OK
PC-50M	Type A	OK	OK
PET-508MA	Type B	–	–
PET-512MA/512MC	Type B	–	–
PET-511BTM* ¹	Type C	–	OK
PET-512MD* ¹ , PEI-512VX* ¹	Type C	–	–
PET-609MA* ¹	Type C	–	–
PET-805LA	Type C	–	–
PLT-604AT/704AT/805AT/1204BT PVT-375BT/375SC/382BT/674BT	Type D	OK	OK
PST-25BT/65AT	Type D	OK	OK
PVT-745BTF/745BTH	Type D	–	OK
PLT-308P/308BTP	Type D	OK	OK
PVT-770RT	Type D	–	OK
PLT-705BTF/705BTH	Type D	OK* ²	OK
PVT-745BTV	Type D	OK* ²	OK
PLT-1202S	Type D	OK* ²	OK
PVT-350BTP	Type D	–	OK
PLT-1204BX	Type D	OK	OK
PLT-704SBT/705BT/1005BT	Type D	OK	OK
PVT-781VT/781VTE	Type D	–	OK
PVT-661VT	Type D	–	OK
PLT-704ST	Type D	OK	OK
PLT-1204ST, PVT-375ST	Type D	OK	OK

OK : Use of the chemical is permitted. : Use of the chemical is not permitted.

*1: If a waterproof cover is placed over the transducer connector, the entire transducer can be immersed in water or cleaning solution.

*2: Only if the transducer is not used internally.

Benzalkonium chloride/ Didecyldimethyl- ammonium chloride/ Polyhexamethylene biguanide (PHMB)	Propan-2-ol Ethanol	Didecyldimethylammonium chloride/Alkyl (C12-16) dimethylbenzylammonium chloride/Alkyl (C12-C14) ethylbenzylammonium chloride	Propan-1-ol Ethanol	
Clinell® Universal Wipes	mikrozid® universal wipes premium	mikrozid® sensitive wipes	mikrozid® AF liquid	mikrozid® AF wipes
GAMA Healthcare Ltd.	Schulke & Mayr GmbH	Schulke & Mayr GmbH	Schulke & Mayr GmbH	Schulke & Mayr GmbH
GBR	DEU	DEU	DEU	DEU
Working solution	Working solution	Working solution	Working solution	Working solution
Wiping only	Wiping only	Wiping only	Wiping only	Wiping only
Room temperature	Room temperature	Room temperature	Room temperature	Room temperature
Normal humidity	Normal humidity	Normal humidity	Normal humidity	Normal humidity
Normal pressure	Normal pressure	Normal pressure	Normal pressure	Normal pressure
-	-	-	-	-
-	-	-	-	-
-	-	-	-	-
-	-	-	-	-
-	-	-	-	-
-	-	-	-	-
-	-	-	-	-
-	-	-	-	-
-	-	-	-	-
OK	OK	OK	OK	OK
OK	OK	OK	OK	OK
-	-	OK	OK	OK
-	-	OK	OK	OK
-	-	-	-	-
-	-	-	-	-
-	-	-	-	-
-	-	OK	OK	OK
-	-	-	-	-
OK	OK	OK	OK	OK
OK	OK	OK	OK	OK
OK	OK	OK	OK	OK
-	-	-	-	-
-	-	-	-	-
-	-	-	-	-

Low/Middle-Level Disinfection -1 (2/2) (Except for USA)

Note that local regulations may require that the transducer be disinfected before sterilization.

Country of origin :
 Any : Any country
 FRA : France
 DEU : Germany
 GBR : United Kingdom
 USA : United States
 JPN : Japan
 AUS : Australia
 CAN : Canada

	Chemical name/type	Didecyldimethyl-ammonium chlorides/Quaternary ammonium compounds	Didecyldimethyl-ammonium chlorides	Glyoxal/Glutaral/Didecyldimethyl-ammonium chlorides
	Trade name	CLEANISEPT® WIPES	WIP'ANIOS Excel	Lysoformin® 3000
	Manufacturer	Dr. Schumacher GmbH	Laboratories ANIOS	Lysoform Disinfection Ltd.
	Country of origin	DEU	FRA	DEU
	Concentration (Dilution ratio)	Working solution	Working solution	2% (50)
	Time	Wiping only	Wiping only	15 min.
	Temperature	Room temperature	Room temperature	Room temperature
	Humidity	Normal humidity	Normal humidity	Normal humidity
	Pressure	Normal pressure	Normal pressure	Normal pressure
Model name	Pressure	Normal pressure	Normal pressure	Normal pressure
PST-25SX	Immersible range	-	-	-
	Type D			
PST-50BT, PVT-712BT	Type D	OK	OK	OK
PST-25ST	Type D	OK	-	OK
PST-30BT	Type D	OK	OK	OK
PLT-1204MV PVT-375MV/382MV/675MV/675MVL	Type D	OK	OK	OK
PVT-681MVL	Type D	-	OK	OK
PVT-681MV	Type D	-	OK	-
PVL-715RST	Type D	-	OK	-
PVI-475BX/574BX/482BX PVT-574BT, PLI-705BX PSI-30BX, PST-28BT/65BT	Type D	OK	-	-
PLI-1205BX/2004BX/605BX PVI-475BT, PVT-475BT/482BT	Type D	OK	-	-
PSI-70BT	Type D	OK	OK	OK
PSI-30VX/40VX	Type D	OK	-	-
PSI-50VX	Type D	OK	-	-
PLI-2002BT	Type D	OK	-	-
PVI-450BXP	Type D	OK	-	-
PLI-3003BX	Type D	OK	-	OK
PLT-1202BT	Type D	OK	OK	OK
PVT-675MVS	Type D	OK	OK	-

OK : Use of the chemical is permitted. : Use of the chemical is not permitted.

*1: If a waterproof cover is placed over the transducer connector, the entire transducer can be immersed in water or cleaning solution.

*2: Not for use in the USA.

Benzalkonium chloride/ Didecyldimethyl- ammonium chloride/ Polyhexamethylene biguanide (PHMB)	Propan-2-ol Ethanol	Didecyldimethylammonium chloride/Alkyl (C12-16) dimethylbenzylammonium chloride/Alkyl (C12-C14) ethylbenzylammonium chloride	Propan-1-ol Ethanol	
Clinell® Universal Wipes	mikrozid® universal wipes premium	mikrozid® sensitive wipes	mikrozid® AF liquid	mikrozid® AF wipes
GAMA Healthcare Ltd.	Schulke & Mayr GmbH	Schulke & Mayr GmbH	Schulke & Mayr GmbH	Schulke & Mayr GmbH
GBR	DEU	DEU	DEU	DEU
Working solution	Working solution	Working solution	Working solution	Working solution
Wiping only	Wiping only	Wiping only	Wiping only	Wiping only
Room temperature	Room temperature	Room temperature	Room temperature	Room temperature
Normal humidity	Normal humidity	Normal humidity	Normal humidity	Normal humidity
Normal pressure	Normal pressure	Normal pressure	Normal pressure	Normal pressure
-	-	-	-	-
OK	OK	OK	OK	OK
-	-	-	-	-
OK	OK	OK	OK	OK
-	-	-	-	-
OK	OK	OK	OK	OK
-	-	-	-	-
OK	OK	OK	OK	OK
OK	OK	OK	OK	OK
OK	-	OK	-	-
OK	OK	OK	OK	OK
OK	OK	OK	OK	OK
OK	OK	OK	OK	OK
-	-	OK	-	-
-	-	OK	OK	OK
-	-	-	-	-
-	-	OK	OK	OK
OK	-	-	-	-

Low/Middle-Level Disinfection -2 (1/2)

Note that local regulations may require that the transducer be disinfected before sterilization.

Country of origin :
 Any : Any country
 FRA : France
 DEU : Germany
 GBR : United Kingdom
 USA : United States
 JPN : Japan
 AUS : Australia
 CAN : Canada

	Chemical name/type	Ethanol	Isopropyl alcohol	Isopropyl alcohol/ Ethylene glycol monobutyl ether	Sodium hypochlorite
Trade name	Ethanol	Isopropyl alcohol	CaviWipes™	Dispatch® Hospital Cleaner Disinfectant Towels with Bleach	
Manufacturer	Any manufacturer	Any manufacturer	Metrex Research	The Clorox Company	
Country of origin	Any	Any	USA	USA	
Concentration (Dilution ratio)	80% (original solution)	70% (original solution)	Working solution	0.65% (working solution)	
Time	Wiping only	Wiping only	Wiping only	Wiping only	
Temperature	Room temperature	Room temperature	Room temperature	Room temperature	
Humidity	Normal humidity	Normal humidity	Normal humidity	Normal humidity	
Model name	Pressure	Normal pressure	Normal pressure	Normal pressure	Normal pressure
PC-20M	Immersible range	–	OK	OK	OK
	Type A	–	OK	OK	OK
PC-50M	Type A	OK	OK	–	–
PET-508MA	Type B	OK	–	–	–
PET-512MA/512MC	Type B	OK*2	–	–	–
PET-511BTM*1	Type C	OK	–	–	–
PET-512MD*1, PEI-512VX*1	Type C	OK	OK	–	–
PET-609MA*1	Type C	OK	–	–	–
PET-805LA	Type C	–	–	–	–
PLT-604AT/704AT/805AT/1202S/1204BT PVT-375BT/375SC/382BT/674BT	Type D	OK	OK	OK	OK
PST-25BT/65AT	Type D	OK	OK	OK	OK
PVT-745BTF/745BTH	Type D	OK	OK	OK	OK
PLT-308P/308BTP	Type D	OK	OK	OK	OK
PVT-770RT	Type D	OK	OK	OK	OK
PLT-705BTF/705BTH	Type D	OK	OK	OK	–
PVT-745BTV	Type D	OK	OK	OK	–
PVT-350BTP	Type D	OK	OK	–	–
PLT-1204BX	Type D	OK	OK	OK	OK
PLT-704SBT/705BT/1005BT PVT-781VT/781VTE	Type D	OK	OK	OK	OK
PLT-704ST, PVT-661VT	Type D	OK	OK	OK	OK
PLT-1204ST, PVT-375ST	Type D	OK	OK	OK	OK

OK : Use of the chemical is permitted. : Use of the chemical is not permitted.

*1: If a waterproof cover is placed over the transducer connector, the entire transducer can be immersed in water or cleaning solution.

*2: Not for use in the USA.

Quaternary ammonium chlorides		Quaternary ammonium chlorides/ Isopropyl alcohol		Decyldimethyloctyl- ammonium chloride/ Dimethyldioctyl- ammonium chloride/ Didecyldimethyl- ammonium chloride/ Quaternary ammonium compounds, benzyl-C12-16- alkyldimethyl, chlorides	Hydrogen peroxide	Hydrogen peroxide	Isopropanol/ Ethanol/ Ethylene glycol monobutyl ether (2- Butoxyethanol)/ Didecyldimethyl- ammonium chloride/ Water	Octyl decyl dimethyl ammonium chloride
Protex™ Ultra Disinfectant Wipes	Sani-Cloth® AF3, Sani-Cloth® AF	Sani-Cloth® Plus	Super Sani-Cloth®	SONO Ultrasound Wipes	Oxivir 1 wipes	PREempt Wipes	CaviWipes1™	Trophon® Companion Cleaning Wipes
Parker Laboratories, Inc.	Professional Disposables International, Inc.	Professional Disposables International, Inc.		Advanced Ultrasound Solutions Inc.	Diversey, Inc.	Virox Technologies Inc.	Metrex Research	Nanosonics Limited
USA	USA	USA		USA	USA	CAN	USA	AUS
0.232% (working solution)	0.28% (working solution)	0.25% / 14.85% (working solution)	0.5% / 55% (working solution)	Working solution	0.5% (working solution)	0.5% (working solution)	Working solution	Working solution
Wiping only	Wiping only	Wiping only	Wiping only	Wiping only	Wiping only	Wiping only	Wiping only	Wiping only
Room temperature	Room temperature	Room temperature	Room temperature	Room temperature	Room temperature	Room temperature	Room temperature	Room temperature
Normal humidity	Normal humidity	Normal humidity	Normal humidity	Normal humidity	Normal humidity	Normal humidity	Normal humidity	Normal humidity
Normal pressure	Normal pressure	Normal pressure	Normal pressure	Normal pressure	Normal pressure	Normal pressure	Normal pressure	Normal pressure
OK	OK	OK	OK	-	OK	OK	-	-
-	OK	OK	OK	-	-	-	-	-
-	OK	-	-	-	-	-	-	-
-	OK	-	-	-	-	-	-	-
-	-	-	-	-	-	-	-	-
-	OK	-	-	-	-	-	-	-
-	OK	-	-	-	-	-	-	-
-	-	-	-	-	-	-	-	-
OK	OK	OK	OK	OK	OK	OK	OK	OK
OK	OK	OK	OK	OK	OK	OK	OK	OK
OK	OK	OK	OK	OK	OK	OK	OK	OK
OK	OK	OK	OK	OK	OK	OK	OK	OK
OK	OK	OK	OK	-	OK	OK	-	-
-	OK	OK	OK	-	-	-	-	-
-	OK	OK	OK	-	-	-	-	-
-	-	-	-	-	-	-	-	-
OK	OK	OK	OK	OK	OK	OK	OK	OK
OK	OK	OK	OK	OK	OK	OK	OK	OK
OK	OK	OK	OK	-	OK	-	-	-
OK	OK	OK	OK	-	OK	-	-	-

Low/Middle-Level Disinfection -2 (2/2)

Note that local regulations may require that the transducer be disinfected before sterilization.

Country of origin :
 Any : Any country
 FRA : France
 DEU : Germany
 GBR : United Kingdom
 USA : United States
 JPN : Japan
 AUS : Australia
 CAN : Canada

	Chemical name/type	Ethanol	Isopropyl alcohol	Isopropyl alcohol/ Ethylene glycol monobutyl ether	Sodium hypochlorite
Trade name	Ethanol	Isopropyl alcohol	CaviWipes™	Dispatch® Hospital Cleaner Disinfectant Towels with Bleach	
Manufacturer	Any manufacturer	Any manufacturer	Metrex Research	The Clorox Company	
Country of origin	Any	Any	USA	USA	
Concentration (Dilution ratio)	80% (original solution)	70% (original solution)	Working solution	0.65% (working solution)	
Time	Wiping only	Wiping only	Wiping only	Wiping only	
Temperature	Room temperature	Room temperature	Room temperature	Room temperature	
Humidity	Normal humidity	Normal humidity	Normal humidity	Normal humidity	
Model name	Pressure	Normal pressure	Normal pressure	Normal pressure	Normal pressure
PST-25SX	Immersible range	OK	OK	-	-
	Type D	OK	OK	OK	OK
PST-50BT, PVT-712BT	Type D	OK	OK	OK	OK
PST-25ST	Type D	OK	OK	OK	OK
PST-30BT	Type D	OK	OK	OK	OK
PLT-1204MV, PVT-375MV/382MV PVT-675MV/675MVL	Type D	OK	OK	OK	OK
PVT-681MVL	Type D	OK	OK	OK	OK
PVT-681MV	Type D	OK	OK	OK	OK
PVL-715RST	Type D	OK	OK	OK	OK
PVI-475BX/475BT/574BX/482BX PVT-475BT/482BT/574BT PLI-1205BX/2004BX/705BX/605BX PSI-30BX, PST-28BT/65BT	Type D	OK	OK	OK	OK
PVI-450BXP	Type D	OK	OK	OK	OK
PSI-70BT	Type D	OK	OK	OK	OK
PSI-30VX/40VX	Type D	OK	OK	OK	OK
PSI-50VX	Type D	OK	OK	OK	OK
PLI-2002BT	Type D	OK	OK	OK	OK
PLI-3003BX	Type D	OK	OK	OK	-
PLT-1202BT	Type D	OK	OK	OK	OK
PVT-675MVS	Type D	OK	OK	OK	OK

OK : Use of the chemical is permitted. : Use of the chemical is not permitted.

*1: If a waterproof cover is placed over the transducer connector, the entire transducer can be immersed in water or cleaning solution.

*2: Not for use in the USA.

High-Level Disinfection -1 (1/2) (Except for USA)

Note that local regulations may require that the transducer be disinfected before sterilization.

Country of origin :

Any : Any country

FRA : France

DEU : Germany

GBR : United Kingdom

USA : United States

JPN : Japan

AUS : Australia

CAN : Canada

	Chemical name/type	Hydrogen peroxide, Peracetic acid	Chlorine dioxide			
	Trade name	PeraSafe™ *2	Tristel Fuse for Instruments	Tristel Sporicidal Wipes/ Tristel Rinse Wipes	Tristel Duo ULT	
	Manufacturer	DuPont	Tristel Solutions Limited			
	Country of origin	Any	GBR			
	Concentration (Dilution ratio)	1.62% (61.7)	0.012% (1 sachet in 5 liters of water)	Working solution		
	Time	10 min.	5 min.	Wiping only		
	Temperature	Room temperature	Room temperature			
	Humidity	Normal humidity	Normal humidity			
	Model name	Pressure	Normal pressure	Normal pressure		
		Immersible range				
PC-20M	Type A	OK	OK	OK	OK	
PC-50M	Type A	OK	OK	OK	OK	
PET-508MA	Type B	OK	–	OK	–	
PET-512MA/512MC	Type B	OK	–	OK	–	
PET-511BTM*1	Type C	OK	OK	OK	–	
PET-512MD*1, PEI-512VX*1	Type C	OK	–	OK	–	
PET-609MA*1	Type C	OK	–	OK	–	
PET-805LA	Type C	–	–	–	–	
PLT-604AT/704AT/805AT/1202S/1204BT PVT-375BT/375SC/382BT/661VT/674BT	Type D	OK	OK	OK	OK	
PST-25BT/65AT	Type D	OK	OK	OK	OK	
PVT-745BTF/745BTH	Type D	OK	OK	OK	OK	
PLT-308P/308BTP, PVT-770RT	Type D	OK	OK	OK	OK	
PLT-705BTF/705BTH	Type D	OK	OK	OK	OK	
PVT-745BTV	Type D	OK	OK	OK	OK	
PVT-350BTP	Type D	OK	OK	–	–	
PLT-1204BX	Type D	OK	OK	OK	OK	
PLT-704SBT/705BT/1005BT PVT-781VT/781VTE	Type D	OK	OK	OK	OK	
PLT-704ST	Type D	OK	OK	OK	OK	
PLT-1204ST, PVT-375ST	Type D	OK	OK	OK	OK	

OK : Use of the chemical is permitted. : Use of the chemical is not permitted.

*1: If a waterproof cover is placed over the transducer connector, the entire transducer can be immersed in water or cleaning solution.

*2: PeraSafe™ may be classified as a sterilizing agent or as a high-level disinfectant depending on the local regulations.

*3: For bacteria and fungi

*4: For virucidal

*5: For mycobactericidal

Succindialdehyde		Peracetic acid			Ortho-phthalaldehyde	Sodium percarbonate/ Citric acid	Peracetic acid, Hydrogen peroxide, Acetic acid			
gigasept® FF (new)		NU-CIDEX®	Anioxyde 1000		Opal™	Clinell® Sporicidal wipes	gigasept® PAA concentrate		mikrozyd® PAA wipes	
Schülke & Mayr GmbH		J&J	Laboratories ANIOS		Whiteley Medical	GAMA Healthcare Ltd.	Schülke & Mayr GmbH			
DEU		GBR	FRA		AUS	GBR	DEU			
5% (20)		0.35% (working solution)	0.15% (working solution)		0.57% (Working solution)	Working solution	2% (a double-chambered bottle into 5 liters of water)		Working solution	
15 min.*3	60 min.*4	5 min.	30 min.	5 min.	6 min.	Wiping only		5 min.*5	15 min.*4	Wiping only
Room temperature		Room temperature			Room temperature	Room temperature	Room temperature			
Normal humidity		Normal humidity			Normal humidity	Normal humidity	Normal humidity			
Normal pressure		Normal pressure			Normal pressure	Normal pressure	Normal pressure			
OK	-	OK	OK	-	-	-	OK	-	-	
-	-	-	OK	-	OK	-	OK	-	-	
OK	OK	-	-	-	OK	OK	OK	OK	OK	
OK	OK	-	-	-	OK	OK	OK	OK	OK	
OK	OK	OK	OK	-	OK	-	OK	OK	OK	
OK	OK	-	-	-	OK	OK	OK	OK	OK	
OK	OK	OK	-	-	OK	OK	OK	OK	OK	
-	-	-	-	OK	OK	-	-	-	-	
OK	OK	OK	OK	-	OK	-	OK	OK	OK	
OK	OK	OK	OK	-	OK	-	OK	OK	OK	
OK	OK	-	OK	-	OK	-	OK	OK	OK	
OK	OK	OK	OK	-	OK	-	OK	OK	OK	
OK	OK	-	OK	-	OK	-	OK	-	OK	
OK	OK	-	OK	-	OK	-	OK	-	OK	
OK	-	-	OK	-	OK	-	OK	-	-	
OK	OK	-	OK	-	OK	-	OK	OK	OK	
OK	OK	-	OK	-	OK	-	OK	OK	OK	
OK	OK	OK	OK	-	OK	-	OK	OK	OK	
OK	OK	-	OK	-	OK	-	OK	OK	OK	
OK	OK	OK	OK	-	OK	-	OK	OK	OK	

High-Level Disinfection - 1 (2/2) (Except for USA)

Note that local regulations may require that the transducer be disinfected before sterilization.

Country of origin :

Any : Any country

FRA : France

DEU : Germany

GBR : United Kingdom

USA : United States

JPN : Japan

AUS : Australia

CAN : Canada

Model name	Chemical name/type	Hydrogen peroxide, Peracetic acid	Chlorine dioxide		
	Trade name	PeraSafe™ *2	Tristel Fuse for Instruments	Tristel Sporicidal Wipes/ Tristel Rinse Wipes	Tristel Duo ULT
	Manufacturer	DuPont	Tristel Solutions Limited		
	Country of origin	Any	GBR		
	Concentration (Dilution ratio)	1.62% (61.7)	0.012% (1 sachet in 5 liters of water)	Working solution	
	Time	10 min.	5 min.	Wiping only	
	Temperature	Room temperature	Room temperature		
	Humidity	Normal humidity	Normal humidity		
	Pressure	Normal pressure	Normal pressure		
	Immersible range	OK	–	–	–
PST-25SX	Type D	OK	–	–	–
PST-25ST/50BT, PVT-712BT	Type D	OK	OK	OK	OK
PST-30BT	Type D	OK	OK	OK	OK
PLT-1204MV PVT-375MV/382MV/675MV/675MVL	Type D	OK	OK	–	–
PVT-681MV	Type D	–	–	–	–
PVT-681MVL	Type D	OK	OK	OK	OK
PVL-715RST	Type D	OK	OK	OK	OK
PVI-475BX/475BT/574BX/482BX/450BXP PVT-475BT/482BT/574BT PLI-1205BX/2004BX/705BX/605BX PSI-30BX, PST-28BT/65BT	Type D	OK	OK	OK	OK
PSI-70BT	Type D	OK	OK	OK	OK
PSI-30VX/40VX	Type D	OK	OK	OK	OK
PSI-50VX	Type D	OK	OK	OK	OK
PLI-2002BT	Type D	OK	OK	OK	OK
PLI-3003BX	Type D	–	OK	–	OK
PLT-1202BT	Type D	OK	OK	OK	OK
PVT-675MVS	Type D	OK	–	OK	OK

OK : Use of the chemical is permitted. : Use of the chemical is not permitted.

*1: If a waterproof cover is placed over the transducer connector, the entire transducer can be immersed in water or cleaning solution.

*2: PeraSafe™ may be classified as a sterilizing agent or as a high-level disinfectant depending on the local regulations.

*3: For bacteria and fungi

*4: For virucidal

*5: For mycobactericidal

Succindialdehyde		Peracetic acid		Ortho-phthalaldehyde	Sodium percarbonate/ Citric acid	Peracetic acid, Hydrogen peroxide, Acetic acid		
gigasept® FF (new)		NU-CIDEX®	Anioxyde 1000	Opal™	Clinell® Sporidical wipes	gigasept® PAA concentrate		mikrozyd® PAA wipes
Schülke & Mayr GmbH		J&J	Laboratories ANIOS	Whiteley Medical	GAMA Healthcare Ltd.	Schülke & Mayr GmbH		
DEU		GBR	FRA	AUS	GBR	DEU		
5% (20)		0.35% (working solution)	0.15% (working solution)	0.57% (Working solution)	Working solution	2% (a double-chambered bottle into 5 liters of water)		Working solution
15 min.*3	60 min.*4	5 min.	30 min.	6 min.	Wiping only	5 min.*5	15 min.*4	Wiping only
Room temperature		Room temperature		Room temperature	Room temperature	Room temperature		
Normal humidity		Normal humidity		Normal humidity	Normal humidity	Normal humidity		
Normal pressure		Normal pressure		Normal pressure	Normal pressure	Normal pressure		
OK	-	-	-	OK	-	OK	-	-
OK	OK	OK	OK	OK	-	OK	OK	OK
OK	OK	OK	OK	OK	-	OK	OK	OK
OK	-	OK	OK	OK	-	OK	-	OK
OK	OK	-	-	OK	-	OK	OK	OK
OK	OK	OK	OK	OK	-	OK	OK	OK
OK	OK	-	OK	OK	-	OK	OK	OK
-	-	OK	OK	OK	-	OK	OK	OK
OK	OK	OK	OK	OK	-	OK	OK	OK
-	-	OK	OK	OK	-	OK	OK	OK
-	-	OK	OK	OK	-	OK	OK	OK
-	-	OK	OK	OK	-	OK	OK	OK
-	-	OK	OK	OK	-	-	-	-
OK	OK	OK	OK	OK	-	OK	OK	OK
-	-	OK	OK	OK	-	-	-	OK

High-Level Disinfection -2 (1/2)

Note that local regulations may require that the transducer be disinfected before sterilization.

Country of origin :
 Any : Any country
 FRA : France
 DEU : Germany
 GBR : United Kingdom
 USA : United States
 JPN : Japan
 AUS : Australia
 CAN : Canada

Chemical name/type	Glutaraldehyde				
	Trade name	CIDEX [®]	CIDEX PLUS [®] 28 day solution	MetriCide [™]	MetriCide [™] 28
Manufacturer	J&J		Metrex Research		
Country of origin	Any		USA		
Concentration (Dilution ratio)	2.4% (working solution)	3.4% (working solution)	2.6% (working solution)	2.5% (working solution)	3.4% (working solution)
Time	45 min.	20 min.	45 min.	90 min.	90 min.
Temperature	Room temperature		25°C (77°F)		
Humidity	Normal humidity		Normal humidity		
Pressure	Normal pressure		Normal pressure		
Model name	Pressure				
PC-20M	Immersible range	OK		OK	
	Type A	OK		OK	
PC-50M	Type A	OK		-	
PET-508MA	Type B	OK		OK	
PET-512MA/512MC	Type B	OK		OK	
PET-511BTM ^{*1}	Type C	OK		OK	
PET-512MD ^{*1} , PEI-512VX ^{*1}	Type C	OK	OK ^{*3}	OK	
PET-609MA ^{*1}	Type C	OK	OK ^{*3}	OK	
PET-805LA	Type C	-		-	
PLT-604AT/704AT/805AT/1202S/1204BT PVT-375BT/375SC/382BT/674BT	Type D	OK		OK	
PST-25BT/65AT	Type D	OK		OK	
PVT-745BTF/745BTH	Type D	OK		OK	
PLT-308P/308BTP	Type D	OK		OK	
PVT-770RT	Type D	OK		OK	
PLT-705BTF/705BTH	Type D	OK		OK	
PVT-745BTV	Type D	OK		OK	
PVT-350BTP	Type D	OK		OK	
PLT-1204BX	Type D	OK		OK	
PLT-704SBT/705BT/1005BT PVT-781VT/781VTE	Type D	OK		OK	
PLT-704ST, PVT-661VT	Type D	OK		OK	
PLT-1204ST, PVT-375ST	Type D	OK		OK	

OK : Use of the chemical is permitted. : Use of the chemical is not permitted.

*1: If a waterproof cover is placed over the transducer connector, the entire transducer can be immersed in water or disinfectant.

*2: For transducers with a thin cable, use a Trophon[®] EPR with a serial number beginning with 24100 or higher or a transducer cable clamp "T-Clamp". Please contact your Trophon[®] EPR representative for further information.

*3: Not for use in the USA.

*4: It is recommended to use Trophon[®] Companion Wipes or mikrozid[®] sensitive wipes for the reprocessing requirements of the ultrasound transducers, including when Trophon[®] EPR/Trophon[®] 2 is used for high-level disinfection.

Glutaraldehyde		Ortho-phthalaldehyde	Hydrogen peroxide			
WAVICIDE® -01	Sporicidin®	CIDEX® OPA	Trophon® EPR*4	Trophon® 2*4	SPOROX® II	Revital-Ox™ Resert® High Level Disinfectant
Medical Chemical Corporation	Contec Inc.	J&J	Nanosonics Limited		Sultan Healthcare	STERIS Corporation
USA	USA	Any	AUS		USA	USA
2.5% (working solution)	1.12% (working solution)	0.55% (working solution)	Use the disinfectant cartridges specifically designed for the disinfection device. There are no disinfection related parameters (temperature, humidity, pressure, or time) to be set by the operator.		7.5% (working solution)	2% (working solution)
45 min.	20 min.	12 min.			30 min.	8 min.
Room temperature	Room temperature	Room temperature			Room temperature	Room temperature
Normal humidity	Normal humidity	Normal humidity			Normal humidity	Normal humidity
Normal pressure	Normal pressure	Normal pressure			Normal pressure	Normal pressure
OK	-	-	-	-	OK	OK
OK	-	OK	-	-	-	-
OK	-	OK	-	-	-	-
OK	-	OK	-	-	-	-
OK	OK	OK	-	-	-	-
OK	-	OK	-	-	-	OK
OK	-	OK	-	-	-	OK
-	-	OK	-	-	-	-
OK	OK	OK	OK	OK	OK	OK
OK	OK	OK	OK*2	OK	OK	OK
OK	OK	OK	OK*2	OK	OK	OK
OK	OK	OK	-	-	OK	OK
OK	OK	OK	OK	OK	OK	OK
OK	OK	OK	OK	OK	OK	-
OK	OK	OK	OK*2	OK	OK	-
OK	OK	OK	-	-	-	-
OK	OK	OK	OK	OK	OK	OK
OK	OK	OK	OK	OK	OK	OK
OK	OK	OK	OK	OK	OK	OK
OK	OK	OK	OK	-	OK	OK
OK	OK	OK	OK	-	OK	OK

High-Level Disinfection -2 (2/2)

Note that local regulations may require that the transducer be disinfected before sterilization.

Country of origin :

Any : Any country

FRA : France

DEU : Germany

GBR : United Kingdom

USA : United States

JPN : Japan

AUS : Australia

CAN : Canada

Chemical name/type	Glutaraldehyde				
	Trade name	CIDEX PLUS [®] 28 day solution	MetriCide™	MetriCide™ 28	MetriCide™ Plus 30
Manufacturer	J&J		Metrex Research		
Country of origin	Any		USA		
Concentration (Dilution ratio)	2.4% (working solution)	3.4% (working solution)	2.6% (working solution)	2.5% (working solution)	3.4% (working solution)
Time	45 min.	20 min.	45 min.	90 min.	90 min.
Temperature	Room temperature		25°C (77°F)		
Humidity	Normal humidity		Normal humidity		
Pressure	Normal pressure		Normal pressure		
Model name	Immersible range				
PST-25SX	Type D	OK		-	
PST-50BT, PVT-712BT	Type D	OK		OK	
PST-25ST	Type D	OK		OK	
PST-30BT	Type D	OK		OK	
PLT-1204MV, PVT-375MV/382MV PVT-675MV/675MVL	Type D	OK		OK	
PVT-681MVL	Type D	OK		OK	
PVT-681MV	Type D	OK		-	
PVL-715RST	Type D	OK		OK	
PVI-475BT/574BX, PVT-475BT/482BT/574BT PLI-1205BX/2004BX/605BX	Type D	OK		OK	
PVI-475BX/482BX/450BXP, PSI-30BX PST-28BT/65BT, PLI-705BX	Type D	OK		OK	
PSI-70BT	Type D	OK		OK	
PSI-30VX/40VX	Type D	OK		OK	
PSI-50VX	Type D	OK		OK	
PLI-2002BT	Type D	OK		OK	
PLI-3003BX	Type D	OK		OK	
PLT-1202BT	Type D	OK		OK	
PVT-675MVS	Type D	OK		OK	

OK : Use of the chemical is permitted. : Use of the chemical is not permitted.

- *1: If a waterproof cover is placed over the transducer connector, the entire transducer can be immersed in water or disinfectant.
- *2: For transducers with a thin cable, use a Trophon[®] EPR with a serial number beginning with 24100 or higher or a transducer cable clamp "T-Clamp". Please contact your Trophon[®] EPR representative for further information.
- *3: Not for use in the USA.
- *4: It is recommended to use Trophon[®] Companion Wipes or mikrozid[®] sensitive wipes for the reprocessing requirements of the ultrasound transducers, including when Trophon[®] EPR/Trophon[®] 2 is used for high-level disinfection.

Glutaraldehyde		Ortho-phthalaldehyde	Hydrogen peroxide			
WAVICIDE® -01	Sporicidin®	CIDEX® OPA	Trophon® EPR*4	Trophon® 2*4	SPOROX® II	Revital-Ox™ Resert® High Level Disinfectant
Medical Chemical Corporation	Contec Inc.	J&J	Nanosonics Limited		Sultan Healthcare	STERIS Corporation
USA	USA	Any	AUS		USA	USA
2.5% (working solution)	1.12% (working solution)	0.55% (working solution)	Use the disinfectant cartridges specifically designed for the disinfection device. There are no disinfection related parameters (temperature, humidity, pressure, or time) to be set by the operator.		7.5% (working solution)	2% (working solution)
45 min.	20 min.	12 min.			30 min.	8 min.
Room temperature	Room temperature	Room temperature			Room temperature	Room temperature
Normal humidity	Normal humidity	Normal humidity			Normal humidity	Normal humidity
Normal pressure	Normal pressure	Normal pressure			Normal pressure	Normal pressure
-	OK	OK	-	-	-	-
OK	OK	OK	OK*2	OK	OK	OK
OK	OK	OK	OK*2	-	OK	OK
OK	OK	OK	OK*2	OK	OK	OK
OK	-	OK	OK	OK	OK	OK
OK	-	OK	OK	OK	OK	OK
-	-	OK	OK	OK	OK	OK
OK	OK	OK	-	-	OK	OK
OK	OK	OK	OK*2	OK	OK	OK
OK	OK	OK	OK*2	OK	OK	OK
OK	OK	OK	OK*2	OK	OK	OK
OK	OK	OK	OK*2	OK	OK	OK
OK	OK	OK	OK*2	OK	OK	OK
OK	OK	OK	OK*2	OK	OK	OK
OK	OK	OK	-	-	-	-
OK	OK	OK	OK*2	OK	OK	OK
OK	-	OK	OK	OK	OK	OK

Sterilization

Note that local regulations may require that the transducer be disinfected before sterilization.

Country of origin :
 Any : Any country
 FRA : France
 DEU : Germany
 GBR : United Kingdom
 USA : United States
 JPN : Japan
 AUS : Australia
 CAN : Canada

Model name	Pressure	Chemical name/type		Ethylene oxide gas ^{*2}		Hydrogen peroxide plasma ^{*3}		
		Trade name	Manufacturer	Country of origin	Concentration (Dilution ratio)	Time	Temperature	Humidity
PC-20M/50M	980 hPa [gauge]	Ethylene oxide gas	Any manufacturer	Any	10%	20% ^{*7}	J&J	Any
PET-508MA ^{*1} /512MA ^{*1} /512MC ^{*1} /511BTM ^{*1}							STERRAD [®] 50/100S/200/NX	STERRAD [®] 100NX
PET-512MD ^{*1} /609MA ^{*1} , PEI-512VX ^{*1}							Use STERRAD system cassettes specifically designed for the sterilizer.	
PET-805LA							Select [SHORT Cycle] for STERRAD 100S/200 ^{*5} and [STANDARD Cycle] for STERRAD NX.	Select [EXPRESS Cycle]. There are no sterilization related parameters (temperature, humidity, pressure, or time) to be set by the operator.
PLT-604AT/704AT/805AT, PST-25BT/65AT		OK					There are no sterilization related parameters (temperature, humidity, pressure, or time) to be set by the operator.	There are no sterilization related parameters (temperature, humidity, pressure, or time) to be set by the operator.
PLT-1204BT, PVT-382BT/661VT/674BT		OK						
PVT-375BT/375SC		OK						OK ^{*6, *8}
PVT-745BTF/745BTH		OK				OK ^{*9}	OK	OK ^{*6, *8}
PLT-308P ^{*4} , PVT-770RT		OK					OK	
PLT-705BTF/705BTH		OK				OK ^{*9}	OK	OK ^{*6, *8}
PVT-745BTV		OK				OK ^{*9}	OK	OK ^{*6, *8}
PVT-350BTP ^{*4} , PLT-308BTP ^{*4}		OK				OK ^{*9}	OK	OK ^{*6, *8}
PLT-1204BX		OK					OK	
PLT-704SBT/705BT/1005BT							OK	OK ^{*6, *8}
PLT-704ST								
PLT-1204ST, PVT-375ST								
PST-25SX								
PST-25ST/50BT, PVT-712BT								
PST-30BT		OK						
PLT-1204MV, PVT-375MV/382MV PVT-675MV/675MVL/681MV/681MVL								
PVT-781VT/781VTE						OK ^{*9}	OK	OK ^{*6, *8}
PVL-715RST							OK	
PVI-475BX/475BT/574BX/482BX PVT-475BT/482BT/574BT PLI-1205BX/2004BX/705BX/605BX PSI-30BX, PST-28BT/65BT						OK ^{*9}	OK	OK ^{*6, *8}
PSI-70BT							OK	OK ^{*6, *8}
PLT-1202S		OK					OK	OK ^{*6, *8}
PSI-30VX/40VX								OK ^{*6, *8}
PSI-50VX								
PLI-2002BT						OK ^{*9}	OK	OK ^{*6, *8}
PLI-3003BX								
PLT-1202BT						OK ^{*9}	OK	OK ^{*6, *8}
PVT-675MVS								
PVI-450BXP						OK ^{*10}	OK	OK ^{*6, *8}


OK : Use of the chemical is permitted. - : Use of the chemical is not permitted.

-
- *1: The mouthpieces supplied with the TEE transducer can be sterilized by high-pressure steam sterilization (temperature: 134°C (273.2°F), holding time: 18 min). Confirm that "AUTOCLAVABLE" is indicated on the mouthpiece.
 - *2: After sterilization, thoroughly degas the transducer to remove all gas residues on the transducer.
 - *3: Note that the label on the transducer may fade and the connector handle may become difficult to turn. These are not abnormalities.
 - *4: The biopsy adaptors, needle stopper, and depth gauge supplied with the transducer can be also sterilized by high-pressure steam sterilization (temperature: 134°C (273.2°F), holding time: 18 min).
 - *5: For STERRAD 100S/200 for the USA, no Cycle option is provided.
 - *6: One transducer can be sterilized per sterilization cycle. Do not perform sterilization for more than one item at a time.
 - *7: Not for use in the USA.
 - *8: There may be sticky areas on the cable, however, this is not a fault. If there are sticky areas, wipe them with a soft cloth or gauze moistened with ethanol.
 - *9: These items have been confirmed to withstand 120 cycles of sterilization using the method described in this manual.
 - *10: These items have been confirmed to withstand 240 cycles of sterilization using the method described in this manual.

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<https://eu.medical.canon>

For other countries, please contact your local distributor.

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