



# OS and 3rd Party Patch List

## X-ray

	FMI Name	Marketing Name	Security Patch Applicable Version	Microsoft OS Security patch identification information
<b>XR</b>	FMI-OSXR-CVE-2019-0708	HDR-08A/P7,/G7,/S7,/B7 HDR-08A/P7 (Ultimax-i) and S7 (Kalare) with V1.60.04 and earlier HDR-08A/P7 with V1.60.04 and later	HDR-08A	KB4499175 for CVE-2019-0708
<b>XR</b>	CXDI Windows 10 LTSB 64bit January 2021	<a href="#">KB4576750</a>	Windows 10 Version 1607 Service Stack Update	
		<a href="#">KB890830</a>	Malicious Software Removal Tool	
		<a href="#">KB4346087</a>	Windows 10 Version 1607 Update	
		<a href="#">KB4580325</a>	Windows 10 Version 1607 Security Update for Adobe Flash Player	
		<a href="#">KB4598243</a>	Cumulative update for Windows 10 Version 1607	
		<a href="#">KB4535680</a>	Security Update for Windows 10 Version 1607	
	"CXDI Windows 10 LTSB 64bit April 2021"	<a href="#">KB5001402</a>	2021-04 Servicing Stack Update for Windows 10 Version 1607 for x64-based Systems	
		<a href="#">KB4346087</a>	2019-02 Update for Windows 10 Version 1607 for x64-based Systems	
		<a href="#">KB4580325</a>	2020-10 Security Update for Adobe Flash Player for Windows 10 Version 1607 for x64-based Systems	
		<a href="#">KB4535680</a>	Security Update for Windows 10 Version 1607 for x64-based Systems	
		<a href="#">KB5001347</a>	2021-04 Cumulative Update for Windows 10 Version 1607 for x64-based Systems	
		<a href="#">KB4577586</a>	Update for Removal of Adobe Flash Player for Windows 10 Version 1607 for x64-based systems	

**CANON MEDICAL SYSTEMS USA, INC.**

2441 Michelle Drive, Tustin CA 92780 | 800.421.1968



	FMI Name	Marketing Name	Security Patch Applicable Version	Microsoft OS Security patch identification information
XR	CXDI Windows 10 LTSCB	<a href="#">KB5001402</a>	2021-04 Servicing Stack Update for Windows 10 Version 1607 for x64-based Systems	
		<a href="#">KB4346087</a>	2019-02 Update for Windows 10 Version 1607 for x64-based Systems	
		<a href="#">KB4535680</a>	Security Update for Windows 10 Version 1607 for x64-based Systems	
		<a href="#">KB4589210</a>	2021-01 Update for Windows 10 Version 1607 for x64-based Systems	
		<a href="#">KB4577586</a>	Update for Removal of Adobe Flash Player for Windows 10 Version 1607 for x64-based systems	
		<a href="#">KB5004238</a>	2021-07 Cumulative Update for Windows 10 Version 1607 for x64-based Systems	
XR	FMI-OSXR-APR2021	Microsoft Windows 10 Enterprise LTSC 2019 (x64)	"HDR-08A/PX (Ultimax-i) V2.00.02 HDR-08A/BX (PlessartZERO) V2.00.02 HDR-08A/SX V2.00.02 HDR-08A/KX V2.00.02"	<a href="#">FMI-OSXR-APR2021.pdf</a>
XR	CXDI Windows10 LTSC 2019 64 bit Oct 2021	<a href="#">KB4598480</a>	2021-01 Servicing Stack Update for Windows Server 2019 for x64-based Systems	
		<a href="#">KB5004424</a>	2021-07 Servicing Stack Update for Windows Server 2019 for x64-based Systems	
		<a href="#">KB4465065</a>	2019-02 Update for Windows Server 2019 for x64-based Systems	
		<a href="#">KB4577586</a>	Update for Removal of Adobe Flash Player for Windows Server 2019 for x64-based systems	
		<a href="#">KB5006765</a>	2021-10 Cumulative Update for .NET Framework 3.5, 4.7.2 and 4.8 for Windows Server 2019 for x64	
		<a href="#">KB4535680</a>	Security Update for Windows Server 2019 for x64-based Systems	
		<a href="#">KB4589208</a>	2021-01 Update for Windows Server 2019 for x64-based Systems	
		<a href="#">KB5006744</a>	2021-10 Cumulative Update Preview for Windows Server 2019 for x64-based Systems	
		<a href="#">KB890830</a>	Windows Malicious Software Removal Tool x64 - v5.94	

**CANON MEDICAL SYSTEMS USA, INC.**

2441 Michelle Drive, Tustin CA 92780 | 800.421.1968



	FMI Name	Marketing Name	Security Patch Applicable Version	Microsoft OS Security patch identification information
<b>XR</b>	CXDI Windows10 LTSB 2016 64 bit Oct 2021	<a href="#">KB5005698</a>	2021-09 Servicing Stack Update for Windows Server 2016 for x64-based Systems	
		<a href="#">KB4346087</a>	2019-02 Update for Windows Server 2016 for x64-based Systems	
		<a href="#">KB4535680</a>	Security Update for Windows Server 2016 for x64-based Systems	
		<a href="#">KB4589210</a>	2021-01 Update for Windows Server 2016 for x64-based Systems	
		<a href="#">KB4577586</a>	Update for Removal of Adobe Flash Player for Windows Server 2016 for x64-based systems	
		<a href="#">KB5006669</a>	2021-10 Cumulative Update for Windows Server 2016 for x64-based Systems	
		<a href="#">KB890830</a>	Windows Malicious Software Removal Tool x64 - v5.94	

**CANON MEDICAL SYSTEMS USA, INC.**

2441 Michelle Drive, Tustin CA 92780 | 800.421.1968

Canon Medical Systems USA, INC. 2022. All rights reserved. Design and specifications are subject to change without notice.

*Made For life*