



# OS and 3rd Party Patch List

## X-ray

	FMI Name	Marketing Name	Security Patch Applicable Version	Microsoft OS Security patch identification information
<b>XR</b>	FMI-OSXR-CVE-2019-0708	HDR-08A/P7,/G7,/S7,/B7 HDR-08A/P7 (Ultimax-i) and S7 (Kalare) with V1.60.04 and earlier HDR-08A/P7 with V1.60.04 and later	HDR-08A	KB4499175 for CVE-2019-0708
<b>XR</b>	CXDI Windows 10 LTSB 64bit January 2021	<a href="#">KB4576750</a>	Windows 10 Version 1607 Service Stack Update	
		<a href="#">KB890830</a>	Malicious Software Removal Tool	
		<a href="#">KB4346087</a>	Windows 10 Version 1607 Update	
		<a href="#">KB4580325</a>	Windows 10 Version 1607 Security Update for Adobe Flash Player	
		<a href="#">KB4598243</a>	Cumulative update for Windows 10 Version 1607	
		<a href="#">KB4535680</a>	Security Update for Windows 10 Version 1607	
	"CXDI Windows 10 LTSB 64bit April 2021"	<a href="#">KB5001402</a>	2021-04 Servicing Stack Update for Windows 10 Version 1607 for x64-based Systems	
		<a href="#">KB4346087</a>	2019-02 Update for Windows 10 Version 1607 for x64-based Systems	
		<a href="#">KB4580325</a>	2020-10 Security Update for Adobe Flash Player for Windows 10 Version 1607 for x64-based Systems	
		<a href="#">KB4535680</a>	Security Update for Windows 10 Version 1607 for x64-based Systems	
		<a href="#">KB5001347</a>	2021-04 Cumulative Update for Windows 10 Version 1607 for x64-based Systems	
		<a href="#">KB4577586</a>	Update for Removal of Adobe Flash Player for Windows 10 Version 1607 for x64-based systems	

**CANON MEDICAL SYSTEMS USA, INC.**

2441 Michelle Drive, Tustin CA 92780 | 800.421.1968



	FMI Name	Marketing Name	Security Patch Applicable Version	Microsoft OS Security patch identification information
XR	CXDI Windows 10 LTSCB	<a href="#">KB5001402</a>	2021-04 Servicing Stack Update for Windows 10 Version 1607 for x64-based Systems	
		<a href="#">KB4346087</a>	2019-02 Update for Windows 10 Version 1607 for x64-based Systems	
		<a href="#">KB4535680</a>	Security Update for Windows 10 Version 1607 for x64-based Systems	
		<a href="#">KB4589210</a>	2021-01 Update for Windows 10 Version 1607 for x64-based Systems	
		<a href="#">KB4577586</a>	Update for Removal of Adobe Flash Player for Windows 10 Version 1607 for x64-based systems	
		<a href="#">KB5004238</a>	2021-07 Cumulative Update for Windows 10 Version 1607 for x64-based Systems	
XR	FMI-OSXR-APR2021	Microsoft Windows 10 Enterprise LTSC 2019 (x64)	"HDR-08A/PX (Ultimax-i) V2.00.02 HDR-08A/BX (PlessartZERO) V2.00.02 HDR-08A/SX V2.00.02 HDR-08A/KX V2.00.02"	<a href="#">FMI-OSXR-APR2021.pdf</a>
XR	CXDI Windows10 LTSC 2019 64 bit Oct 2021	<a href="#">KB4598480</a>	2021-01 Servicing Stack Update for Windows Server 2019 for x64-based Systems	
		<a href="#">KB5004424</a>	2021-07 Servicing Stack Update for Windows Server 2019 for x64-based Systems	
		<a href="#">KB4465065</a>	2019-02 Update for Windows Server 2019 for x64-based Systems	
		<a href="#">KB4577586</a>	Update for Removal of Adobe Flash Player for Windows Server 2019 for x64-based systems	
		<a href="#">KB5006765</a>	2021-10 Cumulative Update for .NET Framework 3.5, 4.7.2 and 4.8 for Windows Server 2019 for x64	
		<a href="#">KB4535680</a>	Security Update for Windows Server 2019 for x64-based Systems	
		<a href="#">KB4589208</a>	2021-01 Update for Windows Server 2019 for x64-based Systems	
		<a href="#">KB5006744</a>	2021-10 Cumulative Update Preview for Windows Server 2019 for x64-based Systems	
		<a href="#">KB890830</a>	Windows Malicious Software Removal Tool x64 - v5.94	

**CANON MEDICAL SYSTEMS USA, INC.**

2441 Michelle Drive, Tustin CA 92780 | 800.421.1968



	FMI Name	Marketing Name	Security Patch Applicable Version	Microsoft OS Security patch identification information	
<b>XR</b>	CXDI Windows10 LTSB 2016 64 bit Oct 2021	<a href="#">KB5005698</a>	2021-09 Servicing Stack Update for Windows Server 2016 for x64-based Systems		
		<a href="#">KB4346087</a>	2019-02 Update for Windows Server 2016 for x64-based Systems		
		<a href="#">KB4535680</a>	Security Update for Windows Server 2016 for x64-based Systems		
		<a href="#">KB4589210</a>	2021-01 Update for Windows Server 2016 for x64-based Systems		
		<a href="#">KB4577586</a>	Update for Removal of Adobe Flash Player for Windows Server 2016 for x64-based systems		
		<a href="#">KB5006669</a>	2021-10 Cumulative Update for Windows Server 2016 for x64-based Systems		
			<a href="#">KB890830</a>	Windows Malicious Software Removal Tool x64 - v5.94	
	CXDI Windows10 LTSC 2019 64 bit Aug 2022	<a href="#">KB4465065</a>	2019-02 Update for Windows 10 Version 1809 for x64-based Systems		
		<a href="#">KB4598480</a>	2021-01 Servicing Stack Update for Windows 10 Version 1809 for x64-based Systems		
		<a href="#">KB4535680</a>	Security Update for Windows 10 Version 1809 for x64-based Systems		
		<a href="#">KB5005112</a>	2021-08 Servicing Stack Update for Windows 10 Version 1809 for x64-based Systems		
		<a href="#">KB5012647</a>	2022-04 Cumulative Update for Windows 10 Version 1809 for x64-based Systems		
		<a href="#">KB5012128</a>	2022-04 Cumulative Update for .NET Framework 3.5 and 4.7.2 for Windows 10 Version 1809 for x64		
			<a href="#">KB890830</a>	Windows Malicious Software Removal Tool x64 - v5.99	
	CXDI Windows10 LTSB 2016 64 bit Aug 2022	<a href="#">KB5011570</a>	2022-03 Servicing Stack Update for Windows 10 Version 1607 for x64-based Systems		
		<a href="#">KB4535680</a>	Security Update for Windows 10 Version 1607 for x64-based Systems		
		<a href="#">KB4589210</a>	2021-01 Update for Windows 10 Version 1607 for x64-based Systems		
		<a href="#">KB4577586</a>	Update for Removal of Adobe Flash Player for Windows 10 Version 1607 for x64-based systems		
<a href="#">KB5012596</a>		2022-04 Cumulative Update for Windows 10 Version 1607 for x64-based Systems			
		<a href="#">KB890830</a>	Windows Malicious Software Removal Tool x64 - v5.99		

**CANON MEDICAL SYSTEMS USA, INC.**

2441 Michelle Drive, Tustin CA 92780 | 800.421.1968



	FMI Name	Marketing Name	Security Patch Applicable Version	Microsoft OS Security patch identification information		
<b>XR</b>	CXDI Windows10 V21H2 64 bit Jan 2023	<a href="#">KB890830</a>	Windows Malicious Software Removal Tool x64 - v5.109			
		<a href="#">KB5020872</a>	2022-12 Cumulative Update for .NET Framework 3.5 and 4.8 for Windows 10 Version 21H2 for x64			
		<a href="#">KB5022282</a>	2023-01 Cumulative Update for Windows 10 Version 21H2 for x64-based Systems			
		<a href="#">KB5012170</a>	2022-08 Security Update for Windows 10 Version 21H2 for x64-based Systems			
	CXDI Windows10 LTSC 2019 64bit Jan 2023	<a href="#">KB890830</a>	Windows Malicious Software Removal Tool x64 - v5.109			
		<a href="#">KB4052623</a>	Update for Microsoft Defender Antivirus antimalware platform - KB4052623 (Version 4.18.2211.5)		<a href="https://catalog.s.download.windowsupdate.com/d/msdownload/update/software/defu/2022/11/updateplat-form_b5a2679b058450feb68b78736e525f8f5ac657fb.exe">https://catalog.s.download.windowsupdate.com/d/msdownload/update/software/defu/2022/11/updateplat-form_b5a2679b058450feb68b78736e525f8f5ac657fb.exe</a>	
					<a href="https://catalog.s.download.windowsupdate.com/d/msdownload/update/software/defu/2022/11/updateplat-form_038cfdedde8500f926be3462b3e41a9be3e76436.exe">https://catalog.s.download.windowsupdate.com/d/msdownload/update/software/defu/2022/11/updateplat-form_038cfdedde8500f926be3462b3e41a9be3e76436.exe</a>	
		<a href="#">KB4465065</a>	2019-02 Update for Windows 10 Version 1809 for x64-based Systems			
		<a href="#">KB4577586</a>	Update for Removal of Adobe Flash Player for Windows 10 Version 1809 for x64-based systems			
		<a href="#">KB5009472</a>	2022-02 Cumulative Update for .NET Framework 3.5 and 4.7.2 for Windows 10 Version 1809 for x64			
		<a href="#">KB5017270</a>	2022-10 Cumulative Update for .NET Framework 3.5 and 4.7.2 for Windows 10 Version 1809 for x64			
		<a href="#">KB5020866</a>	2022-12 Cumulative Update for .NET Framework 3.5 and 4.7.2 for Windows 10 Version 1809 for x64			
		<a href="#">KB5005112</a>	2021-08 Servicing Stack Update for Windows 10 Version 1809 for x64-based Systems			
		<a href="#">KB4589208</a>	2021-01 Update for Windows 10 Version 1809 for x64-based Systems			
		<a href="#">KB5022286</a>	2023-01 Cumulative Update for Windows 10 Version 1809 for x64-based Systems			
		<a href="#">KB5012170</a>	2022-08 Security Update for Windows 10 Version 1809 for x64-based Systems			
		CXDI Windows10 LTSB2016 64bit Jan 2023	<a href="#">KB4103720</a>	2018-05 Cumulative Update for Windows 10 Version 1607 for x64-based Systems		
			<a href="#">KB5017396</a>	2022-09 Servicing Stack Update for Windows 10 Version 1607 for x64-based Systems		
	<a href="#">KB890830</a>		Windows Malicious Software Removal Tool x64 - v5.109			
	<a href="#">KB4346087</a>		2019-02 Update for Windows 10 Version 1607 for x64-based Systems			
	<a href="#">KB4589210</a>		2021-01 Update for Windows 10 Version 1607 for x64-based Systems			
	<a href="#">KB4577586</a>		Update for Removal of Adobe Flash Player for Windows 10 Version 1607 for x64-based systems			
	<a href="#">KB5012170</a>		2022-08 Security Update for Windows 10 Version 1607 for x64-based Systems			
	<a href="#">KB5022289</a>		2023-01 Cumulative Update for Windows 10 Version 1607 for x64-based Systems			
	<a href="#">KB4052623</a>		Update for Windows Defender Antivirus antimalware platform - (Version 4.18.2001.10)		<a href="https://catalog.s.download.windowsupdate.com/c/msdownload/update/software/updt/2020/02/updateplat-form_0456e6719c3ee098af03b785230ac020643fa1ac.exe">https://catalog.s.download.windowsupdate.com/c/msdownload/update/software/updt/2020/02/updateplat-form_0456e6719c3ee098af03b785230ac020643fa1ac.exe</a>	

**CANON MEDICAL SYSTEMS USA, INC.**

2441 Michelle Drive, Tustin CA 92780 | 800.421.1968



	FMI Name	Marketing Name	Microsoft OS Security patch identification information
XR	CXDI Windows10 V22H2 64 bit Feb 2023	<a href="#">KB890830</a>	Windows Malicious Software Removal Tool x64 - v5.110
		<a href="#">KB4052623</a>	Update for Microsoft Defender Antivirus antimalware platform
		<a href="https://catalog.s.download.windowsupdate.com/c/msdownload/update/software/defu/2023/01/updateplatform_ec690aca321b2b1828f710ba74bdd934764249ad.exe">https://catalog.s.download.windowsupdate.com/c/msdownload/update/software/defu/2023/01/updateplatform_ec690aca321b2b1828f710ba74bdd934764249ad.exe</a>	
		<a href="https://catalog.s.download.windowsupdate.com/c/msdownload/update/software/defu/2023/01/updateplatform_72f81c826b05ab8c3d118c21ae6bfbe7921aec76.exe">https://catalog.s.download.windowsupdate.com/c/msdownload/update/software/defu/2023/01/updateplatform_72f81c826b05ab8c3d118c21ae6bfbe7921aec76.exe</a>	
		<a href="#">KB5012170</a>	2022-12 Security Update for Windows 10 Version 22H2 for x64-based Systems
		<a href="#">KB5022502</a>	2023-02 Cumulative Update for .NET Framework 3.5 and 4.8 for Windows 10 Version 22H2 for x64
		<a href="#">KB5022834</a>	2023-02 Cumulative Update for Windows 10 Version 22H2 for x64-based Systems
	CXDI Windows10 LTSC 2019 64bit Feb 2023	<a href="#">KB890830</a>	Windows Malicious Software Removal Tool x64 - v5.110
		<a href="#">KB4052623</a>	Update for Microsoft Defender Antivirus antimalware platform
		<a href="https://catalog.s.download.windowsupdate.com/c/msdownload/update/software/defu/2023/01/updateplatform_ec690aca321b2b1828f710ba74bdd934764249ad.exe">https://catalog.s.download.windowsupdate.com/c/msdownload/update/software/defu/2023/01/updateplatform_ec690aca321b2b1828f710ba74bdd934764249ad.exe</a>	
		<a href="https://catalog.s.download.windowsupdate.com/c/msdownload/update/software/defu/2023/01/updateplatform_72f81c826b05ab8c3d118c21ae6bfbe7921aec76.exe">https://catalog.s.download.windowsupdate.com/c/msdownload/update/software/defu/2023/01/updateplatform_72f81c826b05ab8c3d118c21ae6bfbe7921aec76.exe</a>	
		<a href="#">KB4465065</a>	2019-02 Update for Windows 10 Version 1809 for x64-based Systems
		<a href="#">KB4577586</a>	Update for Removal of Adobe Flash Player for Windows 10 Version 1809 for x64-based systems
		<a href="#">KB5009472</a>	2022-01 Cumulative Update Preview for .NET Framework 3.5 and 4.7.2 for Windows 10 Version 1809 for x64
		<a href="#">KB5022511</a>	2023-02 Cumulative Update for .NET Framework 3.5 and 4.7.2 for Windows 10 Version 1809 for x64
		<a href="#">KB5005112</a>	2021-08 Servicing Stack Update for Windows 10 Version 1809 for x64-based Systems
<a href="#">KB4589208</a>		2021-01 Update for Windows 10 Version 1809 for x64-based Systems	
<a href="#">KB5022840</a>	2023-02 Cumulative Update for Windows 10 Version 1809 for x64-based Systems		
<a href="#">KB5012170</a>	2022-08 Security Update for Windows 10 Version 1809 for x64-based Systems		

**CANON MEDICAL SYSTEMS USA, INC.**

2441 Michelle Drive, Tustin CA 92780 | 800.421.1968



	FMI Name	Marketing Name	Microsoft OS Security patch identification information
XR	CXDI Windows10 LTSB2016 64bit Feb 2023	<a href="#">KB4103720</a>	2018-05 Cumulative Update for Windows 10 Version 1607 for x64-based Systems
		<a href="#">KB5017396</a>	2022-09 Servicing Stack Update for Windows 10 Version 1607 for x64-based Systems
		<a href="#">KB890830</a>	Windows Malicious Software Removal Tool x64 - v5.110
		<a href="#">KB4346087</a>	2019-02 Update for Windows 10 Version 1607 for x64-based Systems
		<a href="#">KB4589210</a>	2021-01 Update for Windows 10 Version 1607 for x64-based Systems
		<a href="#">KB4577586</a>	Update for Removal of Adobe Flash Player for Windows 10 Version 1607 for x64-based systems
		<a href="#">KB5012170</a>	2022-08 Security Update for Windows 10 Version 1607 for x64-based Systems
	CXDI Windows10 LTSC 2019 64bit Feb 2023	<a href="#">KB5022838</a>	2023-02 Cumulative Update for Windows 10 Version 1607 for x64-based Systems
		<a href="#">KB4052623</a>	Update for Windows Defender Antivirus antimalware platform - *patch for 64bit CPU architecture
		<a href="https://catalog.s.download.windowsupdate.com/c/msdownload/update/software/defu/2023/01/updateplatform_ec690aca321b2b-1828f710ba74bd-d934764249ad.exe">https://catalog.s.download.windowsupdate.com/c/msdownload/update/software/defu/2023/01/updateplatform_ec690aca321b2b-1828f710ba74bd-d934764249ad.exe</a>	*Patch for ARM 64 Architecture
	<a href="https://catalog.s.download.windowsupdate.com/c/msdownload/update/software/defu/2023/01/updateplatform_72f-81c826b05ab8c-3d118c21ae6bfbe-7921aec76.exe">https://catalog.s.download.windowsupdate.com/c/msdownload/update/software/defu/2023/01/updateplatform_72f-81c826b05ab8c-3d118c21ae6bfbe-7921aec76.exe</a>	*Patch for 32bit CPU Architecture	

**CANON MEDICAL SYSTEMS USA, INC.**

2441 Michelle Drive, Tustin CA 92780 | 800.421.1968



	FMI Name	Marketing Name	Microsoft OS Security patch identification information
XR	CXDI Windows11 V22H2 64 bit DEC 2023	<a href="#">KB890830</a>	Windows Malicious Software Removal Tool x64 - v5.119
		<a href="#">KB5007651</a>	Update for Windows Security platform antimalware platform - KB5007651 (Version 1.0.2306.10002)
		<a href="https://catalog.s.download.windowsupdate.com/c/msdownload/update/software/defu/2023/06/securityhealthsetup_0b3ffdaa5a5c5b7be1fab0fabfd383859a046189.exe">https://catalog.s.download.windowsupdate.com/c/msdownload/update/software/defu/2023/06/securityhealthsetup_0b3ffdaa5a5c5b7be1fab0fabfd383859a046189.exe</a>	
		<a href="#">KB5012170</a>	2022-08 Security Update for Windows 11 22H2 for x64-based Systems
		<a href="#">KB5032007</a>	2023-11 Cumulative Update for .NET Framework 3.5 and 4.8.1 for Windows 11, version 22H2 for x64
		<a href="#">KB5033375</a>	2023-12 Cumulative Update for Windows 11 Version 22H2 for x64-based Systems
	CXDI Windows10 LTSC 2019 64bit DEC 2023	<a href="#">KB890830</a>	Windows Malicious Software Removal Tool x64 - v5.119
		<a href="#">KB4465065</a>	2019-02 Update for Windows 10 Version 1809 for x64-based Systems
		<a href="#">KB4577586</a>	Update for Removal of Adobe Flash Player for Windows 10 Version 1809 for x64-based systems
		<a href="#">KB5031984</a>	2023-11 Cumulative Update for .NET Framework 3.5 and 4.7.2 for Windows 10 Version 1809 for x64
		<a href="#">KB5005112</a>	2021-08 Servicing Stack Update for Windows 10 Version 1809 for x64-based Systems
		<a href="#">KB4589208</a>	2021-01 Update for Windows 10 Version 1809 for x64-based Systems
		<a href="#">KB5033371</a>	2023-12 Cumulative Update for Windows 10 Version 1809 for x64-based Systems
		<a href="#">KB4019091</a>	Security Update for SQL Server 2014 Service Pack 1 GDR
		<a href="https://catalog.s.download.windowsupdate.com/c/msdownload/update/software/secu/2017/08/sqls-erver2014-kb4019091-x64_e265e2ee82009cf658d187d27928edb3c824c200.exe">https://catalog.s.download.windowsupdate.com/c/msdownload/update/software/secu/2017/08/sqls-erver2014-kb4019091-x64_e265e2ee82009cf658d187d27928edb3c824c200.exe</a>	
		<a href="#">KB2538243</a>	Security Update for Microsoft Visual C++ 2008 Service Pack 1 Redistributable Package
		<a href="https://catalog.s.download.windowsupdate.com/msdownload/update/software/secu/2011/05/vcred-ist_ia64_a01885c05d23b23831b80675a25026e7ad466ce5.exe">https://catalog.s.download.windowsupdate.com/msdownload/update/software/secu/2011/05/vcred-ist_ia64_a01885c05d23b23831b80675a25026e7ad466ce5.exe</a>	
		<a href="https://catalog.s.download.windowsupdate.com/msdownload/update/software/secu/2011/05/vcred-ist_x64_a7c83077b8a28d409e36316d2d7321fa0ccdb7e8.exe">https://catalog.s.download.windowsupdate.com/msdownload/update/software/secu/2011/05/vcred-ist_x64_a7c83077b8a28d409e36316d2d7321fa0ccdb7e8.exe</a>	
		<a href="#">KB2565063</a>	Security Update for Microsoft Visual C++ 2010 Service Pack 1 Redistributable Package
		<a href="https://catalog.s.download.windowsupdate.com/msdownload/update/software/secu/2011/07/vcred-ist_x86_28c54491be70c38c97849c3d8cfbffd0d3c515cb.exe">https://catalog.s.download.windowsupdate.com/msdownload/update/software/secu/2011/07/vcred-ist_x86_28c54491be70c38c97849c3d8cfbffd0d3c515cb.exe</a>	
<a href="https://catalog.s.download.windowsupdate.com/msdownload/update/software/secu/2011/07/vcred-ist_ia64_4b355f7fb318b4b862a20eed32da7f036fadaeff.exe">https://catalog.s.download.windowsupdate.com/msdownload/update/software/secu/2011/07/vcred-ist_ia64_4b355f7fb318b4b862a20eed32da7f036fadaeff.exe</a>			

**CANON MEDICAL SYSTEMS USA, INC.**

2441 Michelle Drive, Tustin CA 92780 | 800.421.1968



	FMI Name	Marketing Name	Microsoft OS Security patch identification information
XR	CXDI Windows10 LTSB2016 64bit DEC 2023	<a href="#">KB4103720</a>	2018-05 Cumulative Update for Windows 10 Version 1607 for x64-based Systems
		<a href="#">KB890830</a>	Windows Malicious Software Removal Tool x64 - v5.119
		<a href="#">KB5032391</a>	2023-11 Servicing Stack Update for Windows 10 Version 1607 for x64-based Systems
		<a href="#">KB4346087</a>	2019-02 Update for Windows 10 Version 1607 for x64-based Systems
		<a href="#">KB4589210</a>	2021-01 Update for Windows 10 Version 1607 for x64-based Systems
		<a href="#">KB4577586</a>	Update for Removal of Adobe Flash Player for Windows 10 Version 1607 for x64-based systems
		<a href="#">KB5012170</a>	2022-08 Security Update for Windows 10 Version 1607 for x64-based Systems
		<a href="#">KB5033373</a>	2023-12 Cumulative Update for Windows 10 Version 1607 for x64-based Systems
		<a href="#">KB4019091</a>	Security Update for SQL Server 2014 Service Pack 1 GDR
		<a href="https://catalog.s.download.windowsupdate.com/c/msdownload/update/software/secu/2017/08/sqls-erver2014-kb4019091-x64_e265e2ee82009cf658d187d27928edb3c824c200.exe">https://catalog.s.download.windowsupdate.com/c/msdownload/update/software/secu/2017/08/sqls-erver2014-kb4019091-x64_e265e2ee82009cf658d187d27928edb3c824c200.exe</a>	
		<a href="#">KB2538243</a>	Security Update for Microsoft Visual C++ 2008 Service Pack 1 Redistributable Package
		<a href="https://catalog.s.download.windowsupdate.com/msdownload/update/software/secu/2011/05/vcred-ist_ia64_a01885c05d23b23831b80675a25026e7ad466ce5.exe">https://catalog.s.download.windowsupdate.com/msdownload/update/software/secu/2011/05/vcred-ist_ia64_a01885c05d23b23831b80675a25026e7ad466ce5.exe</a>	
		<a href="https://catalog.s.download.windowsupdate.com/msdownload/update/software/secu/2011/05/vcred-ist_x86_470640aa4bb7db8e69196b5edb0010933569e98d.exe">https://catalog.s.download.windowsupdate.com/msdownload/update/software/secu/2011/05/vcred-ist_x86_470640aa4bb7db8e69196b5edb0010933569e98d.exe</a>	
		<a href="#">KB2565063</a>	Security Update for Microsoft Visual C++ 2010 Service Pack 1 Redistributable Package
		<a href="https://catalog.s.download.windowsupdate.com/msdownload/update/software/secu/2011/07/vcred-ist_x86_28c54491be70c38c97849c3d8cfbfd0d3c515cb.exe">https://catalog.s.download.windowsupdate.com/msdownload/update/software/secu/2011/07/vcred-ist_x86_28c54491be70c38c97849c3d8cfbfd0d3c515cb.exe</a>	
<a href="https://catalog.s.download.windowsupdate.com/msdownload/update/software/secu/2011/07/vcred-ist_ia64_4b355f7fb318b4b862a20eed32da7f036fadaeff.exe">https://catalog.s.download.windowsupdate.com/msdownload/update/software/secu/2011/07/vcred-ist_ia64_4b355f7fb318b4b862a20eed32da7f036fadaeff.exe</a>			
FMI-OSXR-APR2025_ UDR_W10LTSC	Microsoft Windows 10 IoT Enterprise LTSC 2021 (x64) operating system.	" <a href="#">UDR-10A/C2 V2.6 SP0000*</a> "	
FMI-OSXR-JUN2025_ UDR_W10LTSC	Windows 10 Enterprise LTSC 2019 (x64)	" <a href="#">UDR-10A V2.2 SP0002*-SP004* V2.2 SP0014* V2.2 SP0013*</a> "	
FMI-OSXR-JUL2025_ HDR	Microsoft Windows 10 Enterprise LTSC 2019 (x64)	" <a href="#">HDR-08A/PX, /SX, /KX, /BX V2.00.01 or later HDR-08A/CX V1.70.00 or later HDR-08A/DX V1.80.00 or later</a> "	
FMI-OSXR-JUL2025_ XRDK	Microsoft Windows 10 Enterprise LTSC 2019 (x64)	" <a href="#">XDRK-EREF V2.00*R001 or later</a> "	

**CANON MEDICAL SYSTEMS USA, INC.**

2441 Michelle Drive, Tustin CA 92780 | 800.421.1968



	FMI Name	OS	Security Patch Applicable Version	System Name	CMSC Patch Number
<b>XR</b>	FMI-OSXR-JUN2026	<a href="#">Microsoft Windows 10 Enterprise LTSC 2019 (x64)</a>	"HDR-08A/PX, /SX, /KX, /BX V2.00.01 or later HDR-08A/CX V1.70.00 or later HDR-08A/DX V1.80.00 or later"	"HDR-08A/PX HDR-08A/SX HDR-08A/KX HDR-08A/BX HDR-08A/CX HDR-08A/DX"	<a href="#">FMI-OSXR-JUN2026</a>

**CANON MEDICAL SYSTEMS USA, INC.**

2441 Michelle Drive, Tustin CA 92780 | 800.421.1968

Canon Medical Systems USA, INC. 2025. All rights reserved. Design and specifications are subject to change without notice.