

ENCORE VeloCT Console Upgrade



**"I need upgrades
that meet my needs for
the future with low dose,
accelerated workflow
and compliance."**

Toshiba's new VeloCT console includes exciting upgrades for your existing Aquilion™ CT.


ADAPTIVE ITERATIVE DOSE REDUCTION 3D (AIDR 3D)

- AIDR 3D is fully integrated into scanning protocols and adaptive to your clinical and image quality needs
- Dose Check software includes dose notification and alert features
- DICOM, IHE structured dose report
- XR-29 compliant

ACCELERATED WORKFLOW

- Advanced hardware for faster image reconstruction speeds with AIDR 3D
- coneXact™ supports double slice reconstruction with no dose penalty for improved 3D and multiplanar images

Improve Clinical and Operational Capabilities



Toshiba gives you a voice. What's yours?



TOSHIBA AMERICA MEDICAL SYSTEMS, INC.

2441 Michelle Drive, Tustin CA 92780 / 800.421.1968 medical.toshiba.com

©Toshiba Medical Systems Corporation 2016. All rights reserved.

Aquilion and coneXact are the trademarks of Toshiba Medical Systems Corporation.

Design and specifications are subject to change without notice.

CTCSS11796US
MCACT0121EB