

## ***Ultrasound***



“Developing rapid imaging scan setup and fast data acquisition technologies will be of main interest and concern in emergency radiology, which will significantly increase imaging capacity and reduce risk of delaying emergent patient management.”<sup>1</sup>

### **Rapid Emergency Exam Setup**

“In the ER, time can be both our greatest ally and our worst enemy. The first hour after trauma or the beginning of symptoms is considered the golden hour.”<sup>2</sup>



# Imaging when and where you need it

## Emergency Exam Start



- Start emergency exam without entering patient identification
- Ideal for trauma and time-sensitive scanning

## Customizable User Interface



- Customizable console and mode-sensitive touch screen to meet clinical needs
- Assign functions and preferences to the keys of your choice

## Xario g-series — Ultrasound Anywhere, Anytime



- Ultra-fast boot up time for time-sensitive exams
- Smart standby mode allows exams to be initiated almost instantly<sup>3</sup>



- High mobility brings ultrasound quickly and easily to any location
- Independent monitor and control panel swivel to help optimize viewing in awkward scanning situations



- Optional extended battery, allowing optimum emergency scanning<sup>4,5</sup>
- Enhanced efficiency and transport with optional cable-free, wireless peripheral devices

<sup>1</sup>Shen D (2021) Grand Challenges in Radiology. Frontiers in Radiology.20 April 2021 1:615138

<sup>2</sup>[www.altusemergency.com/the-importance-of-emergency-radiology/](https://www.altusemergency.com/the-importance-of-emergency-radiology/)

<sup>3</sup>Up to 5 second start-up time from smart standby mode

<sup>4</sup>2 hours standard battery on Xario g-series

<sup>5</sup>Up to 8 hours optional extended battery on Xario 200G and up to 4 hours optional on Xario 100G. 30 minutes optional extended battery on Aplio a-series and Aplio i-series

Follow us: <https://us.medical.canon>



@CanonMedicalUS



Canon Medical Systems USA, Inc.



+CanonMedicalUS

CANON MEDICAL SYSTEMS USA, INC.

<https://us.medical.canon> | 2441 Michelle Drive, Tustin CA 92780 | 800.421.1968

©Canon Medical Systems, USA 2022. All rights reserved. Design and specifications subject to change without notice.

Made for Life is a trademark of Canon Medical Systems Corporation.

ULCSS14269US

*Made For life*