

Canon

RADREX™-i

Digital Radiography Systems



Enhance productivity while you improve patient care.

Providing you with high-quality images, RADREX™-i Digital Radiography Systems include unique, fully integrated capabilities and robust ergonomic features designed for high-intensity/high-efficiency environments.



RexView™

Extends control room functions to table side



RexPanel

Expands system flexibility to facilitate out-of-bucky exams



RexSpeed

Fully automated solutions help eliminate redundant tasks while allowing on-the-fly changes

REXPROTECT: REDUCE X-RAY EXPOSURE—OPTIMIZE IMAGE QUALITY.

Set a new standard of safety and lower the dose—for even your smallest patients.



High Detector Quantum Efficiency

Reduce dose while achieving high-quality images



Beam Quality Filters

Help reduce soft radiation to the patient's skin



Manual Collimation

Allows for further reduction in collimation to avoid unnecessary radiation exposure



Exposure Indicator

Lighted displays provide exposure levels at a glance

“If we could
automate
our processes,
we could
streamline
our workflow.”



For higher productivity at a nonstop pace, at the heart of every RADREX-i System are RexSpeed auto-features that include Auto-Mapping, Auto-Center Stop, Auto-Collimation, Auto-Tracking, and Auto-Park.

REXVIEW: THE ROBUST FEATURE YOUR DEPARTMENT NEEDS— AND PATIENTS APPRECIATE.

RexView provides unequalled patient safety with tableside review:

Increase your time with the patient to help enhance patient safety

You no longer have to wait in the control room for images to be displayed

Return to the patient immediately after the exposure

Accept or reject the image at table side

RADREX™-i Digital

Radiography Systems from Canon Medical Systems.

Industry-leading RADREX-i meets the clinical efficiency and image quality demands that today's healthcare facilities need. With an array of unique, fully integrated capabilities, seamless IT integration, and ergonomic features, RADREX-i Systems address your healthcare challenges and provide a better standard of care for your patients.

Health Care's Triple Aim

With RADREX-i, Canon Medical Systems leads the quest for better health care:

- 1 Improved individual care**
RexView™ provides control room functions at tableside so technologists are next to the patient throughout the exam
- 2 Better health for populations**
Reject Analysis software makes care safer by evaluating previously rejected images, allowing for process improvements to reduce re-scans and the risk of unnecessary radiation exposure
- 3 Lower total costs of care**
When compared to CR and film, near instantaneous image preview in the control room means faster procedures, increasing throughput and efficiency. Additionally, elimination of film reduces cost

Follow us: <https://us.medical.canon>



@CanonMedicalUS



+CanonMedicalUS



Canon Medical Systems USA, Inc.



+CanonMedicalUS

Canon

CANON MEDICAL SYSTEMS USA, INC.

<https://us.medical.canon>

2441 Michelle Drive, Tustin CA 92780 | 800.421.1968

©Canon Medical Systems, USA 2018. All rights reserved.

Design and specifications subject to change without notice.

RADREX, RexView and Made for Life are trademarks of Canon Medical Systems Corporation. Google+ logo and YouTube logo are trademarks of Google Inc. TWITTER, TWEET, RETWEET and the Twitter logo are trademarks of Twitter, Inc. or its affiliates. LinkedIn, the LinkedIn logo, the IN logo and InMail are registered trademarks or trademarks of LinkedIn Corporation and its affiliates in the United States and/or other countries.

XRBR12864US MCAXR0325EBA

Made For life