

OMNERA[®] 300M

Robust Performance.
Flexible Functionality.

Introducing the OMNERA 300M Ceiling-Mounted Manual Positioning Digital Radiographic System that offers a blend of performance and flexibility, making it an ideal choice for a wide range of imaging needs.

Leverage the Canon Medical Advantage



Enjoy unprecedented image quality.

Benefit from over a quarter of a century of Canon's expertise in digital detector technology, which offers high image resolution and high sensitivity with low dose potential.

Versatility and value.

Flexibility in system configurations that can be tailored to a department's unique workflow and budget considerations.


Robust software functionality comes standard.

Canon Medical's CXDI Control Software NE has a comprehensive set of features as part of the standard offering, eliminating the need to invest in costly options.

- Advanced Edge Enhancement for improved visibility of catheters, fine structures and bone
- Free Image Rotation enables the orientation of images to any angle
- Radiation Dose Structured Reporting
- Intuitive multi-touch display operation including zoom, pan, double-tap and list scrolling
- Windows 10 OS with cybersecurity features: Bit Locker and Device Guard

Operator-Friendly Features

The overhead tube crane with patient-side touchscreen console, premium wall stand, and robust patient table help support a smooth workflow and a comfortable exam experience for patients. Experience peace of mind with OMNERA 300M's solid design, ensuring reliability when you need it most.



Comfortable Exams for Better Outcomes

- Overhead tube crane with intuitive patient-side touchscreen console displays SID values, tube angulation and generator options right where you need them
- Wall stand features side and overhead-mounted handgrips to help support patient comfort and positioning during upright exams
- Image a broad range of patients with the completely flat and smooth 36 in. wide and 800 lb. weight capacity table that lowers to 22 in. from the floor for easy patient loading

CANON MEDICAL SYSTEMS USA, INC.

<https://us.medical.canon> | 2441 Michelle Drive, Tustin CA 92780 | 800.421.1968

©Canon Medical Systems, USA 2024. All rights reserved. Design and specifications are subject to change without notice. Made for Life is a trademark of Canon Medical Systems Corporation.

Disclaimer: Some features presented in this brochure may not be commercially available on all systems shown or may require the purchase of additional options. Please contact your local representative from Canon Medical Systems for details.

XRBR14537US

Made For life